## **PCT**

## WORLD INTELLECTUAL PROPERTY ORGANIZATION International Bureau



### INTERNATIONAL APPLICATION PUBLISHED UNDER THE PATENT COOPERATION TREATY (PCT)

CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, IL, IS, JP, KE, KG, KP, KR, KZ, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, UZ, VN,	(51) International Patent Classification <sup>6</sup> : A61B 5/00	A1	(11) International Publication Number: WO 97/1789: (43) International Publication Date: 22 May 1997 (22.05.97
	<ul> <li>(22) International Filing Date: 15 November 1996 (</li> <li>(30) Priority Data: 08/559,032 15 November 1995 (15.11.9)</li> <li>(71) Applicant: MEDIMAP, INC. [US/US]; 6955 Cast #208, Frederick, MD 21703 (US).</li> <li>(72) Inventors: KAUFMANN, Peter; Mossvagen 56, S Jama (SE). BEVING, Hakan; Linnegatan 33-35, Stockholm (SE). OLSSON, Nils, U.; 6955 Cast #208, Frederick, MD 21703 (US).</li> <li>(74) Agents: BRUHN, David, E. et al.; Dorsey &amp; Whitne Pillsbury Center South, 220 South Sixth Street, Min</li> </ul>	15.11.9  15.11.9  15.11.9  15.11.9  15.11.9	CA, CH, CN, CZ, DE, DK, EE, ES, FI, GB, GE, HU, II IS, JP, KE, KG, KP, KR, KZ, LK, LR, LS, LT, LU, LV MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU SD, SE, SG, SI, SK, TJ, TM, TR, TT, UA, UG, UZ, VN ARIPO patent (KE, LS, MW, SD, SZ, UG), Eurasian paten (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM), European paten (AT, BE, CH, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU MC, NL, PT, SE), OAPI patent (BF, BJ, CF, CG, CI, CM GA, GN, ML, MR, NE, SN, TD, TG).  Published  With international search report.  Before the expiration of the time limit for amending the claims and to be republished in the event of the receipt of amendments.

#### (57) Abstract

A process for analyzing for a variety of physical medical disorders, which includes: (a) clinically testing by a common analysis, patients free of such disorders to obtain characterization thereof as numerical data; (b) scaling a matrix of the data by computer; (c) configuring an artificial neural network prescribed by the number of variables and number of disorders; (d) fitting the neural network in a computer to the numerical data according to adjustable network parameters which include the number of neurons in a hidden layer, the number of hidden layers, the type of layer transfer function, and the connecting weights; (e) thereby training the neural network; (f) withdrawing from a new patient a sample or samples and subjecting same to like clinical testing to obtain new numeric data; (g) scaling the new data; (h) feeding the scaled new data to the trained neural network; and (i) automatically obtaining a diagnosis of the new patient with respect to such disorders.

## FOR THE PURPOSES OF INFORMATION ONLY

Codes used to identify States party to the PCT on the front pages of pamphlets publishing international applications under the PCT.

AM	Armenia	GB	United Kingdom	MW	Malawi
AT	Austria	GE	Georgia	MX	Mexico
AU	Australia	GN	Guinea	NE	Niger
ВВ	Barbados	GR	Greece	NL.	Netherlands
BE	Belgium	HU	Hungary	NO	Norway
BF	Burkina Faso	IE	Ireland	NZ	New Zealand
BG	Bulgaria	IT	Italy	PI.	Poland
BJ	Benin	JP	Japan	PT	Portugal
BR	Brazil	KE	Kenya	RO	Romania
BY	Belarus	KG	Kyrgystan	RU	Russian Federation
CA	Canada	KP	Democratic People's Republic	SD	Sudan
CF	Central African Republic		of Korea	SE	Sweden
CG	Congo	KR	Republic of Korea	SG	Singapore
СН	Switzerland	KZ	Kazakhstan	SI	Slovenia
CI	Côte d'Ivoire	LI	Liechtenstein	SK	Slovakia
CM	Cameroon	LK	Sri Lanka	SN	Senegal
CN	China	LR	Liberia	SZ	Swaziland
CS	Czechoslovakia	LT	Lithuania	TD	Chad
CZ	Czech Republic	LU	Luxembourg	TG	Togo
DE	Germany	LV	Latvia	TJ	Tajikistan
DK	Denmark	MC	Monaco	TT	Trinidad and Tobago
EE	Estonia	MD	Republic of Moldova	UA	Ukraine
ES	Spain	MG	Madagascar	UG	Uganda
FI	Finland	ML	Mali	US	United States of America
FR	France	MN	Mongolia	UZ	Uzbekistan
GA	Gabon	MR	Mauritania	VN	Viet Nam

# SELECTIVE DIFFERENTIATING DIAGNOSTIC PROCESS BASED ON BROAD DATA BASES

#### Brief Description Of The Invention

5

10

15

20

25

30

35

A diagnostic tool is described that determines a patient's health as to a variety of potential physical and medical disorders. It employs multivariant analysis such as artificial neural networks and partial least squares discriminant analysis in conjunction with standard clinical diagnostic procedures.

#### Background To The Invention

A physician seeking data about a patient will invariably subject the patient to one or more diagnostic procedures. Blood assays or urine assays, and the like, are ordered and a technician draws blood or procures urine from the patient. These samples are automatically analyzed in an analyzer that can generate considerable amounts of data about the sample, but only the data requested by the physician is fed back to him/her. The quality of the physician's services to the patient depends to a large extent on the data given to the physician and the physician's use of the data. If the data fails to note a symptom of a disorder, then the physician may miss the proper diagnosis. Moreover, physician bias can result in misuse or misunderstanding of data causing diagnosis errors. 1

These limitations in the practice of medicine have caused industries to contrive aids that assist the physician in dealing with the data that diagnostic tools are capable of generating. However, such efforts have been limited to specific disorders. Though the tools aid in predicting or finding a disorder with respect to a patient, the tools are so limited that again the data utilized is a mere fraction of that available.

In addition, automatic analyzers are under-utilized because the equipment is not programmed to generate maximum data utilization because of a lack of physician interest in

10

15

20

25

30

35

WO 97/17891 PCT/US96/18254

the additional data. Most physicians are incapable of interpreting data other than standard disorder characterizations.

According to Boyd, "Perspectives on the Use of Chemometrics in Laboratory Medicine," Clin. Chem., volume 32, number 9, pages 1726-1733 (1986):

"Measurements and their interpretation are central in the practice of modern medicine. Laboratory medicine as a subspecialty has played an important role in the development of new diagnostic tests and in the generation of an ever-increasing number of measurements. Although modern laboratory instruments produce prodigious amounts of data regarding measurable physiological and chemical qualities, the abundance of these data does not supply specific answers to questions about the physiology or pathophysiology of various organ systems being evaluated. Unless powerful methods can be discovered for converting data to answers, the measurements generated by modern instrumentation may raise more interpretational problems than they solve. this climate of increasing complexity of data, the necessity for developing tools to aid in interpreting laboratory data has become ever more critical."

A number of systems deal with the data that physiological analysis generates. For example, Drastal et al., <u>Journal of Medical Systems</u>, volume 6, number 5, pages 433-445 (1982) describe a so-called EXPERT medical consultation system in which computer models generate rule systems used for diagnosis. Such a system deals with single disorder evaluations and fails to maximally utilize the data generated. In this particular paper, the authors fail to report the success or failure rate for the system.

Another expert approach is described by Kinney et al., <u>Journal of Medical Systems</u>, volume 12, number 5, pages 3319-326 (1988). Kinney et al. designate the particular disorder or data collection for evaluation and sort out the evaluation according to a relative scale. Though Kinney et al. have a broader approach to data utilization, it too suffers from the

Leape, JAMA, volume 272, number 23, pages 1851-1857 (December 21, 1994).

WO 97/17891

5

10

15

20

25

30

35

fact that the physician's choice of data viewing limits data utilization.

PCT/US96/18254

An interesting approach that uses a multivariant interpretation of the patient's historical data to ascertain evolving abnormalities, is described by H. B. Slotnick et al., Clinical Chemistry, volume 36, number 5, pages 748-751 (1990).

Baxt, "Use of an Artificial Neural Network for Data Analysis in Clinical Decision-Making: The Diagnosis of Acute Coronary Occlusion," Neural Computation, 2, pages 480-489 (1990) employs an artificial neural network to evaluate the data based on substantial and complex input variables by the physician. As a result, the input is sufficient to allow the physician to diagnose whether the patient has an acute coronary occlusion, drawing into question the need for the use of the artificial neural network, other than to aid in refining the results of the physician's examination.

Snow et al., "Artificial Neural Networks In The Diagnosis And Prognosis Of Prostate Cancer: A Pilot Study," The Journal Of Urology, volume 152, pages 1923-1926 (November 1994), describe the analysis for a single disorder, i.e., prostate cancer, based on measurements of serum prostate specific antigen (PSA) and biopsy results of the patient that led to improved false positive and negative conclusions.

Stensmo et al., "A Mixture Model System for Medical and Machine Diagnosis," Teauro, et al. (Ed.) Advances in Neural Information Processing Systems, volume 7, Morgan Kaufmann, San Mateo, CA (1995) attempt to deal with missing data from an analysis by estimating the parameters of the missing data by the EM algorithm.

There is a need for a method that provides a physician with correct analysis and data, which ascertains a medical disorder characterizable from such data and reports such to the physician without dependence on the foresight of the physician.

There is a need for a method which interprets analytical data in a way that it singles out abnormalities for evalua-

10

15

20

25

tion which are not obvious from the direct evaluation of data.

In addition to the foregoing needs, the medical profession needs a method that enhances accuracy of diagnosis, especially about conditions considered abnormal to the physician.

Because of the capacity of an automatic analyzer to generate from a patient's sample, data extending beyond that normally selected for determining the existence or non-existence of a disorder, it would be desirable to have a method for using that data for characterizing known or unknown abnormalities.

It is an objective of this invention to satisfy such needs and desires, as well as others, in a facile and effective manner.

#### The Invention

This invention relates to a process for analyzing for a variety of physical medical disorders, which involves:

- diagnostically analyzing groups of reference patients using one or more known diagnostic procedures and determining whether the patients possess one or more physical medical disorders;
- 2. clinical testing by a common analysis, patients free of such disorders and those patients that possess such one or more disorders, which analysis distinguishes any such disorder thereby obtaining the characterization of a collection of numerical data of such one or more disorders;
- 3. scaling the matrix of said numerical data in a computer;
- 4. configuring in a computer, an artificial neural network prescribed by the number of variables and the number of disorders, which artificial neural network is electronic data that possesses an input layer, one or more hidden layers and an output layer;
- 35 5. fitting the artificial neural network electronic data in a computer to the numerical data according to adjustable parameters of the artificial neural network, including one or more of

a) the number of neurons in a hidden layer;

b) the number of hidden layers;

5

10

15

20

25

30

35

- c) the type of transfer functions in the layers; and
- d) the weights connecting the neurons;
- 6. whereby the artificial neural network is trained in the computer to automatically provide a analytic model;
  - 7. withdrawing from a new patient that has not been diagnosed for such disorder, a sample or samples of a kind taken from said reference patients and subjecting said diagnostic sample to said clinical testing to obtain new numerical data;
  - 8. automatically scaling in the computer the new numerical data to the numerical data derived from the reference patients;
  - 9. feeding the scaled new numerical data to the trained artificial neural network and the analytic model thereof, and
  - 10.automatically obtaining a diagnostic answer of the new patient with respect to such disorders.

This invention contemplates the use of all forms of clinical, chemical, biochemical and biological diagnoses, e.g., the act or process of identifying or determining the nature and cause of a disease or injury through evaluation of patient history, examination, and review of laboratory data. This includes the diagnostic analysis of bodily fluids, cells, soft and hard tissues, organs, and the like. processes include immunoassays, cell size analysis, chromatographic analysis, spectrographic analysis, electrochemical analysis, pH analysis, hematological analysis, ion selective potentiometry and near infrared reflectance spectroscopy (NIR), surface spectroscopy including ion-scattering spectrometry (ISS), secondary ion mass spectrometry (SIMS), and auger emission spectrometry (AES), nuclear magnetic imaging techniques, nuclear magnetic resonance spectrometry (NMR), chemical elemental analysis, gas analysis including gas chromatography, antibody analysis, blotting techniques, electrophoresis, amino acid sequencing, DNA analysis, flow cytometric analysis, ion conductive plasmamass spectrometry (ICP-MS) and atomic absorption spectrometry (AAS), and the like.

Any diagnostic tool applied in the process of the invention is commonly applied in respect to steps (1)-(10) of the process. That means that the assay used in step (1) is applied in step (2) and the results characterized in the remaining steps are based on such assay.

5

10

15

20

25

30

35

An element of the process is that one intentionally includes in the analysis the diagnoses of healthy patients as well as the diagnoses of afflicted patients. Such provides the basis for comparisons that deal with disorders that are unknown to the patient. For example, a patient complains to the physician of an unspecified illness, which the physician could attribute to a variety of disorder conditions. A sample, such as a blood sample, is taken from the patient and the physician asks for a total comparison with the data from a reference group that is known to possess disorders that exhibit similar symptoms to that of the patient as well as a reference group that has no disorders (i.e., healthy pa-The eventually obtained data that is given to the physician characterizes all deviations from the data of such healthy patients as well as data that correlates with the data of the reference group patients that possess the disor-However, the results can be presented to the physician in terms of best matching reference group, thereby assigning the particular disorder for the physician. If no matching disorder is found, then the physician in flagged as to any abnormality of the sample. Consequently, the physician can rely on an analysis that is devoid of subjective input coupled with a probability of accuracy. This aids the physician in determining whether to rely on the output of the analysis and allows the physician to couple the analysis with the physician's own diagnosis.

Scaling is the operation of translating data from a variety of different units of measure into a common characterization that one can mathematically employ. Scaling can be done by analog or digital means and the translations can be used by recourse to computation by hand or by computer. However, scaling is best achieved by translating the data on a computer via digital conversions. The preferred method is by

autoscaling ("Z transformation") which employs computer programs to automatically digitally translate the data. The scaling removes measurement units. Any distribution, with a certain mean (" $\mu$ ") and variance (" $\sigma$ 2") can be reduced to the same normal or shape distribution with mean 0 and variance 1.

5

10

15

20

25

30

35

Harrington, Anal. Chem, 1994, 66, 802-807, characterizes artificial neural networks as "computational systems that simulate the microstructure of biological nervous systems." He goes on to state:

"These systems have broad application to chemometrics because they are powerful pattern recognizers. Backpropagation neural networks (BNNs) are prevalent in analytical applications. BNNs model relations between sets of input and output data by minimizing errors of the network outputs during training. Information is stored by the strength of the connections between processing layers."

The first step is configuring in a computer, an artificial neural network prescribed by the number of variables and the number of disorders, which artificial neural network is electronic data that possesses an input layer, one or more hidden layers and an output layer. The number of variables gives the number of input neurons and the number of disorders will give the number of output neurons. Determining the number of hidden layers and the number of neurons in the hidden layers is part of the optimization process. Optimization comprises (i) determining the right number of hidden layers, (ii) the right number of neurons in the hidden layers, and (iii) overall training of the network.

Choosing the right numbers of (i) and (ii) above is done by configuring the network by a pre-chosen number (e.g., 1,1 or 1,2) and training this network, calculating the prediction error between the real and the model values, storing this data and choosing the next higher number for the number (i) and (ii) above. In the first step in choosing, the pre-chosen number may be at least one hidden layer and at least two neurons. This new network is trained, and the prediction error calculated and stored. This process is repeated until a minimum prediction error is obtained.

10

15

20

25

30

35

PCT/US96/18254

Training of the artificial neural network is a standard procedure in a BNN which evolves by sending data to the network, passing it through the network, receiving the model values and comparing it to the real values, i.e., evaluating the prediction error between the real and model values. This information is back-propagated through the network so that the network weights can be refined stepwise to give selected minimal prediction error. The readjustment is not done haphazardly. All of these steps can be controlled and determined via a computer program.

The selection of type of transfer function in the layer is dependent on the nature of the input data (range of variables), the nature of the output data and the type of answer For a hidden layer one may use the hyperbolic tangent sigmoid function, and for the output neuron, a logistic sigmoid function. However, other functions exist and may be employed instead. For example, a linear transfer function, a threshold function, binary hard limiter, bipolar hard limiter, and the like functions, can be employed instead.

The artificial neural network model can be employed in similarly treating the diagnoses of a sample or samples withdrawn from a new patient that has not been diagnosed for such disorder. The same scaling factors can be used so that the computer automatically calculates the diagnoses.

## DETAILED DESCRIPTION OF THE INVENTION

This invention relates to the use of computer analyses using computing and analytical programs located on a computer hard drive (or optical, magnetic and other memory units), plug-in board, expansion box of a computer containing a standard 486/Pentium or Apple Mac microprocessor (or a later generation microprocessor), employing a standard monitor screen. The invention can be practiced on micro, mini and large frame computers using any basic or specialized operating system such as DOS, UNIX, Windows 3.1, 3.11 and Window95, Windows NT, MAC, Apple, and the like, operating systems. The invention can also be placed on the analytical instrument proper utilizing internal computing resources.

10

15

20

25

The first step in the process of the invention involves determining the analytical parameters to which the patients For example, the analyses can relate to will be subjected. any physiological, chemical and/or pathological issue, without limitation. The analyses can relate to any aspect or all aspects of blood analyses that determine heart related issues, infectious disease disorders, circulatory disease disorders, neural disorders, genetically related issues such as inborn errors of metabolism, disorders related to environmental issues, and the like. The analyses can be of cells, all manner of hard and soft tissues, urine, feces, semen and tissue biopsies, and the like. The analyses (e.g. of blood gases) can relate to all manner of physiological gas and fluid constituents. Having said all of that, it is obvious that the extension of the invention to such a variety of analyses will result in an abundance of data that transcends that requested or designated by the physician seeking the analyses.

The first step is the collection of the reference groups used in the analyses. This done by

- 1. subjecting a candidate reference patient to examination to determine whether or not the patient suffers from the disorder under investigation and categorization.
- 2. collecting the data from as many patients as possible with the disorder. Because such candidate reference patient is regarded to be abnormal because the patient possesses the disorder, a broad spectrum of data is procured about that patient for use in assessing trends in the reference patients and the new patients.
- For example, the following table is a calibration set of blood analysis data from a number of reference patients and healthy controls for the purpose of testing the blood of reference patients with diagnosed heart disease, each of whom was schedule for bypass surgery, reference patients suffering from dental amalgam syndrome, reference patients suffering from exposure to organic solvents, acrylic latex paints, and finally, a control group of healthy patients:

  Calibration set:

TABLE 1

Patient code	WBC	RBC	HB	HCT	MCV	МСН	MCHC	PLT	LYM	LYC
A1A3	6.3	4.32	12.6	40.1	92.9	29.1	31.3	136	30.8	1.9
A1B3	9.4	5.4	15.4	46.2	89.2	28.6	32.1	186	17.1	1.6
A1G3	6.1	4.88	14.8	50.3	103.2	30.4	29.4	155	40.1	2.5
A2A3	6.4	5.32	15.6	43.9	82.6	29.4	35.6	206	28.4	1.8
A2D3	8.8	5.95	17.5	49.6	83.4	29.4	35.2	195	20.7	1.8
A2E3	4.8	5.55	15.1	43	77.6	27.3	35.1	135	26.1	1.3
A2G3	6.8	5.5	15.4	43.6	79.2	28.1	35.5	208	38.7	2.6
A2H3	6.9	5.88	16.2	47	80	27.6	34.5	228	24.4	1.7
A2J3	14.9	5.99	16.5	47.2	78.8	27.5	34.9	194	19.9	3
A2K3	6.3	5.27	16.1	44.7	84.8	30.6	36.1	128	31.9	2
A2L3	6.5	5.4	16.4	45.8	84.9	30.4	35.8	147	41.8	2.7
A3B3	5.5	5.78	15.8	45.9	76.5	27.5	35.9	156	26.1	1.5
A3D3	9.8	5.31	15	44.1	83.2	28.3	34	186	31	3
A3E3	7.9	5.3	14.8	42.8	80.7	27.9	34.7	226	46.3	3.6
A3F3	6.1	5.48	16.4	45.7	83.4	29.7	35.7	160	23	1.3
A3G3	10.6	5.07	14.5	39.2	77.4	28.5	36.8	268	26.4	2.9
A3H3	6.2	6.35	18.1	51.8	81.6	28.5	35	176	30.2	1.9
A313	8.9	4.72	14.2	39.9	84.6	30.1	35.6	413	16.8	1.5
A3M3	6.4	5.58	16.1	45.4	81.4	28.8	35.4	214	35.9	2.3
Amlgm1	12.5	4.35	17.4	45.1	103.7	40.1	38.7	193	22.5	2.8
Amlgm2	5.9	4.2	13.4	39.7	94.4	31.9	33.8	135	38	2.3
Amlgm3	8	4.24	13.7	40.6	95.6	32.2	33.7	150	33	2.6
Amlam4	6.7	5.07	19.1	46.5	96.6	37.7	41.2	191	37.1	2.5
Amlgm6	6	4.25	14.4	39.6	93.2	34	36.4	112	35.1	2.1
Amlgm7	7	4.55	15.5	41.9	92.1	34.1	37	149	26.8	1.9
Amigm8	14.2	4.81	17	45.6	94.9	35.4	37.4	169	21.5	3.1
Amigm9	8.3	4.35	15.3	41.7	95.7	35.1	36.6	159	38.6	3.2
Amigm10	8.2	5.44	18	49.3	90.7	33.1	36.5	147	33.3	2.7
Amlgm11	10.6	4.58	16.3	44	96.1	35.7	37.2	166	33.3	3.5
Amlgm13	6.2	4.38	14.3	39.1	89.3	32.6	36.5	169	41.1	2.5
Amlgm14	6.1	4.86	17.1	46.2	95.3	35.2	36.9	146	40.9	2.5
Amigm15	8.8	4.85	17.3	48.8	100.5	35.7	35.5	170	28.4	2.5
Amlgm16	5.8	3.87	13.3	37.7	97.2	34.4	35.4	158	36.4	2.1
Amlgm17	5.6	5.21	16.5	48.8	93.8	31.8	33.9	221	25.6	1.4
Amlgm20	5.6	4.89	15.7	44	89.9	32.1	35.7	163	32.4	1.9
Amlgm23	5.1	4.15	11.7	38.4	92.6	28.1	30.4	201	35.5	1.8
Amigm24	7.2	5.13	13.2	43.8	85.5	25.6	30	240	20.1	1.4
Amlgm25	6.8	4.41	12	38.6	87.6	27.2	31	233	28.6	1.9
OrgSolv73	. 4.7	4.45	16.6	42.4	95.3	37.3	39.2	226	31.5	1.5
OrgSolv74	4.8	4.55	15.5	40.4	88.8	34	38.3	200	40.7	2
OrgSolv98	4	5.1	16.6	49.7	96.7	32.4	33.5	226	27	1.1
OrgSolv99	9.7	5.2	16.4	49.1	94.1	31.3	33.3	460	29	2.8
OrgSolv67	4.5	4.23	14.9	38.1	90.2	35.2	39	169	50.8	2.3
OrgSolv68	4.1	4.86	15.9	41.7	85.9	32.8	38.2	178	34	1.4
OrgSolv84	4.1	5.4	16.5	51.2	94.5	30.4	32.2	260	35.1	1.5
OrgSolv85	7.2	4.8	16.6	48.8	100.7	34.2	34	329	29.1	2.1
OrgSolv94	4.2	4.6	15.7	45	98.5	34.2	34.7	222	25.2	1.1
OrgSolv95	3.8		17.2	51	97.9	33	33.7	289	23.2	0.9
OrgSolv96	5.8		15.7	47.5	95.7	31.4	32.9	276	38.3	2.2
OrgSolv70	5.9		16.8	43	90.7	35.3	38.9	218	32.3	1.9
OrgSolv71	4.1	4,54	16	43	94.8	35.3	37.2	157	25.6	1
OrgSolv72	4.8	4.8	16.4	43.4	90.4	34.2	37.9	204	49.3	2.4
OrgSolv64	6.2	4.15	13.1	37.7	91	31.6	34.7	272	37.6	2.3
OrgSolv65	3.8	4.58	15.1	42.3	92.3	33.1	35.8	232	35.7	1.3
OrgSolv66	4.7	4.4	15.3	40.2	91.4	34.8	38	211	45.3	2.1
OrgSolv91	3.7	5.1	16.8	49.2	96.1	32.7	34.1	223	32.9	1.2
OrgSolv92	7.8		16.6	49.5	97.4	32.7	33.6	353	24.2	1.9
OrgSolv92	4.7		16.9	49.4	94.5	32.2	34.2	151	24.2	1.2
CRAMOCAN	4.7	5.2	10.9	45.4	54.5	34.2	37.2		2.4.2	1.2

0 0 1 00	0.0	4.0	46.0	50	103.6	34.9	33.7	222	44.7	1.6
OrgSolv88	3.6	4.8	16.8 14.9	44.9	103.4	34.3	33.1	256	42.6	2.1
OrgSolv89	4.9	4.3		48.8	96.9	33.1	34.1	334	36.8	1.8
OrgSolv90	4.9	5	16.6			32.7	36.8		32.6	1.0
OrgSolv61	6.2	4.88	15.9	43.3	88.7		35.6		41.7	1.7
OrgSolv62	3.9	4.5	14.6	40.9 41.4	91.1 92.4	32.5 32.9	35.6		28.6	1.1
OrgSolv63	4	4.48	14.8		89.2	33.5	37.6		28.0	1.5
OrgSolv76	5.5	4.79	16.1	42.8		34.6	38.7	241	32.4	
OrgSolv77	6.3	4.83	16.8	43.3	89.3					2
OrgSolv78	4.8	4.65	16	41.7	89	34.5	38.7	220	50.9	2.4
OrgSolv60	3.6	5.28	16	43.6	82.4	30.3	36.8	204	34.8	1.2
OrgSolv79	3.5	4.7	15.7	41.7	88.7	33.4	37.6		40.6	1.4
OrgSolv80	4.3	4.34	15.1	39	89.8	34.8	38.8	133	46.8	2
OrgSolv81	5.7	4.7	15	44.4	94.1	31.8	33.8	226	37.9	2.1
OrgSolv82	7.2	4.9	17.1	46.8	94.7	34.5	36.5	461	30.9	2.2
OrgSolv86	5.4	4.6	14.9	45.3	98	32.3	33	339	39.2	2.1
Acryl10	6.7	4.97	14.7	44	88.5	29.7	33.5	242	31.9	2.1
Acryl12	7.5	4.72	15.3	45.4	96.1	32.4	33.7	307	18.5	1.4
Acryl13	5	4.28	14.7	44.2	103.2	34.3	33.2	217	39.1	2
Acryl1	4.8	4.85	14.5	43.8	90.4	30	33.2	328	40	1.9
Acryl2	6.3	6.92	14.4	42.8	87	29.2	33.6	339	21.3	1.3
Acryl22	5.3	5	15.5	46.1	90.5	30.5	33.6	347	36.6	1.9
Acryl24	5.5	4.1	13.4	40.4	98.4	32.8	33.3	192	25.7	1.4
Acryl25	7.6	4.93	15	45	91.3	30.5	33.4	338	32.8	2.5
Acryl26	4.2	3.9	12.9	39.4	101	33.1	32.8	272	32.4	1.4
Acryl4	5.2	4.72	14.3	41.6	89.6	30.3	33.8	308	27.5	1.4
Acryl5	4.1	4.18	13.5	40.6	97	32.2	33.3	210	40.5	1.7
Acryl6	7.6	4.96	14.8	45.1	91	29.8	32.8	312	30	2.2
Acryl42	5.1	4.36	13.4	40.3	92.5	30.6	33.1	318	30.8	1.5
Acryl8	5.4	4.2	13.1	39	92.8	31.3	33.6	324	27.2	1.5
Acryl9	6.5	4.86	15.4	43.8	90	30.3	33.8	255	24.6	1.6
Acryl14	4.9	5.25	15.6	46.6	88.8	29.7	33.4	248	31.3	1.5
Acryl16	7.2	4.94	16	45.5	92.2	32.3	35.1	294	44.5	3.2
Acryl17	9.6	5.19	15.3	46.4	89.3	29.6	33.1	342	22.4	2.2
Acryl30	5.6	4.17	13	40	96	31.2	32.5	242	31	1.7
Acryl31	6	4.79	14.1	43.3	90.4	29.5	32.6	207	35.5	2.2
Acryl33	5.4	4.67	14.4	43.2	92.5	30.9	33.4	274	30.4	1.6
Acryl38	4.6	4.86	14.3	43.1	88.5	29.6	33.4	246	44.3	2
Acryl39	8.8	5.61	16.4	48.9	87.2	29.2	33.5	226	33.7	2.9
Acryl18	6.2	4.68	14.9	44.7	95.6	31.7	33.2	283	25.1	1.5
Acryl19	5.7	4.31	13.4	39.2	92.3	31.1	33.7	269	20.1	1.1
Acryl20	5.6	4.67	14.3	43.3	92.8	30.6	32.9	293	27.5	1.6
aC1	5.5	4.53	11.5	37.8	83.6	25.4	30.4	180	35.7	2
aC2	5.9	4.38	12.5	40.7	92.9	28.5	30.7	175	29.9	1.7
aC3	3.9	4.59	12.8	41.4	90.2	27.9	30.9	153	41.7	1.6
aC4	6.1	5.04	14.4	46	91.3	28.6	31.3	175	30.7	1.9
aC5	6.3	4.94	14.1	44.9	90.8	28.5	31.4	178	39.7	2.5
aC7	2.9	4.74	12.3	40.7	85.6	26	30.3	108	45.9	1.3
aC8	5.2	5.27	14.4	45.7	86.8	27.3	31.9	143	34	1.8
C13	4.3	5.09	16.7	43.2	84.9	32.8	38.7		34	1.4
C14	5.4	4.73	17.9	43.9	93	37.9	40.8	127	36.9	2
C16	6.1	4.78	17.2	43.4	94.2	35.9	38.2	153	38.1	2.3
C18	7.4	4.88	18.6	45	92.2	38.1	41.3		30.8	2.3
C20	7.4	5.25	19.2	47	89.4	36.5	40.8	241	26.9	1.9
C20	4.8	4.86	18.3	45.6	93.8	37.8	40.2	144	28.3	1.4
				47.8	87.9	36.1	41.1	193	37.4	2.5
C17	6.7	5.44	19.6		92.5	37.8	40.9	138	23	1.6
C21	6.9	4.82	18.2	44.6	92.5	35.5	40.3	191	38.7	2.8
C22	7.1	5.08	18	44.7			37.4	226		2.5
rC28	5.8	5.58	17.2	45.9	82.1	31.3 27	35.5	206	43 45.8	2.5
rC31	5.7	6.46	17.4	49.3	76.2					
rC41	5.4	5.75	16.2	45	78.5	28.2	35.9	204	30.2	1.6
rC46	4.1	5.77	15.6	44.3	76.7	26.9	35.1	191	43	1.8

rC47	6	5.53	12	36.8	66.5	21.8	32.8	223	17.5	1.1
rC48	4.7	5.2	15.7	44	84.5	30.2	35.8	184	42.6	2
rC51	6.6	5.69	16.2	45.1	79.2	28.4	35.9	165	25.6	1.7
rC52	4	5.25	14.6	41.3	78.8	27.9	35.4	161	39.4	1.6
rC59	5.3	4.98	14.9	43.1	86.6	30	34.7	173	43.2	2.3
rC61	6.3	6.05	16.6	46.6	77.1	27.5	35.6	230	35.3	2.2
rC68	4	5.39	15.6	43.6	80.9	29	35.9	170	44.6	1.8
rC69	4.7	5.59	16.4	45.8	81.9	29.4	35.9	205	32.9	1.6
rC71	6.2	5.74	15.7	45.2	78.7	27.4	34.9	168	44.4	2.7

TABLE 1 (CONT.)

Patient code	RDW	RDS	PCT	MPV	PDW	PDS	MP	MR	MW
A1A3	11.7	9.4	0.09	6.8	8.8	9	4.4	89	147
A1B3	12.3	9.4	0.12	6.6	8.7	8.9	4.6	85	152
A1G3	10.1	10.1	0.11	7.1	8.2	8.5	5.4	92	144
A2A3	10.6	8.1	0.14	7	8.6	8.8	4.7	91	137
A2D3	10.8	8.1	0.13	6.9	8.3	8.6	5	93	145
A2E3	11	8.2	0.09	6.6	8.3	8.7	5	85	141
A2G3	10.9	8.1	0.14	6.8	8.5	8.7	4.6	88	138
A2H3	11.4	8.3	0.16	7	8.4	8.7	4.9	87	133
A2J3	10.9	8.2	0.126	6.5	8.5	8.7	5.2	87	144
A2K3	10.9	8.2	0.08	6.2	8.3	8.7	4.4	93	139
A2L3	10.1	8.1	0.104	7	8.6	8.7	5.5	95	141
A3B3	10.6	8.1	0.112	7.2	8.6	8.9	5.2	85	145
A3D3	11	8.1	0.128	6.9	8.3	8.6	5.8	93	149
A3E3	10.5	8	0.147	6.5	8.3	8.6	4.8	91	152
A3F3	11.2	8.2	0.105	6.6	8.5	8.7	4.2	92	143
A3G3	11.6	8.2	0.181	6.7	8.2	8.6	4.7	85	142
A3H3	11.3	8.1	0.128	7.3	8.4	8.7	5.6	91	154
A313	11.8	8.1	0.269	6.5	8.5	8.8	4.1	95	167
A3M3	10.9	8.1	0.152	7.1	8.5	8.8	5.4	90	145
Amlgm1	11.1	9.2	0.13	6.6	8.6	8.8	5.3	102	176
Amlgm2	10.8	9.3	0.1	7.7	9	9.2	5.8	92	154
Amlgm3	10.8	9.1	0.12	7.7	8.7	8.9	6.5	96	158
Amlgm4	10.7	9.2	0.12	6.5	8.6	8.8	5	90	174
Amlgm6	9.8	9.2	0.08	7.4	9	9.2	6.2	92	170
Amigm7	9.8	9.1	0.11	7.6	8.7	8.8	6.2	92	176
Amlam8	10.5	9.4	0.12	7	8.6	8.8	5.8	92	172
Amlgm9	10.6	9.4	0.11	7.2	8.4	8.6	5.6	92	166
Amlgm10	11.1	9.3	0.11	7.2	8.6	8.8	6.4	88	168
Amigm11	10.1	9.2	0.12	7.1	8.5	8.8	6.2	94	162
Amlgm13	11.8	9.4	0.13	7.4	8.6	8.8	6.2	86	160
Amigm14	10.3	9.2	0.11	7.2	8.8	9	5.2	92	166
Amlgm15	10.7	9.2	0.11	6.7	8.5	8.8	5.2	98	166
Amlgm16	9.8	9.1	0.11	7.1	8.9	9.1	6.4	96	162
Amlgm17	11.2	9.3	0.15	6.8	8.9	9.1	4.6	90	160
Amlgm20	10.6	9.2	0.12	7.2	8.4	8.7	6.4	88	160
Amigm23	10.8	9.2	0.14	7.2	8.4	8.7	5.4	90	152
Amigm24	11.9	9.3	0.17	7.1	8.5	8.7	5.8	84	154
Amlgm25	10.8	9.2	0.17	7.2	7.9	8.3	6	86	150
OrgSolv73	9.9	8.6	0.141	6.2	10.1	10	4.5	99	172
OrgSolv74	10.3	8.7	0.13	6.5	9.8	9.8	5	92	113
OrgSolv98	10.6	9.6	0.158	7	9.8	9.8	5.9	91	153
OrgSolv99	10.5	9.6	0.328	7.1	9.2		5.2	88	151
OrgSolv67	10.3	8.7	0.104	6.2	9.7	9.6	4.8	93	93
OrgSolv68	10.7	10.3	0.125	7.1	10.4	10.3	6.2	89	161
OrgSolv84	10.8	9.5	0.166	6.4	9.2	9.3	5.1	89	149
OrgSolv85	10.8	9.6	0.191	5.8	9.1	9.3	4.6	95	163
OrgSolv94	9.8	9.6	0.148	6.7	9.1	9.2	5.7	92	152
OrgSolv95	10.9		0.188	6.5	9.1	9.2	5	91	155
OrgSolv96	10.4		0.183	6.6	9.2	9.2	5.1	91	137

OrgSolv70	10.6	8.7	0.15	6.9	10.6	10.4	5.5	94	171
OrgSolv71	9.3	8.6		6.6	<del></del>		5		166
OrgSolv72	9.8	8.5	0.118	5.9	9.3		4.5		91
OrgSolv64	10.6	8.6	0.194	7.1			6.2		153
OrgSolv65	10.3	8.7	0.138	5.9			2.9		127
OrgSolv66	9.8	8.6	0.146	6.9		9.7	5		95
OrgSolv91	10.8	9.8		6.6		9.5	5		145
OrgSolv92	10.7	9.5	0.223	6.3	8.7	8.9	5.1	93	149
OrgSolv93	10.7	9.7	0.106	7	10.5		4.5		149
OrgSolv88	10.3	9.5	0.158	7.1	9.7	9.7	5.7	98	111
OrgSolv89	9.8	9.4	0.186	7.3	9.2	9.2	6.7	99	87
OrgSolv90	10.1	9.5	0.213	6.4	8.8	9	5.1	92	132
OrgSolv61	10.5	8.9	0.148	6.4	10.1	10	4.3	90	149
OrgSolv62	10.1	8.8	0.155	6.5	10.2	10.1	5.2	93	113
OrgSolv63	10.1	8.6	0.154	6.7	10.1	10	6	97	143
OrgSolv76	10.3	8.7	0.087	7.6	10.5	11	2.7	92	162
OrgSolv77	10	8.7	0.152	6.3	9.1	9.2	5.3	93	173
OrgSolv78	9.7	8.6	0.137	6.2	9.3	9.4	5	93	91
OrgSolv60	11	8.9	0.129	6.4	10	10	5.2	83	158
OrgSolv79	10.2	8.6	0.075	6.8	10.5	10.3	5.4	93	133
OrgSolv80	10.4	8.8	0.098	7.3	10.5	10.4	6.5	92	89
OrgSolv81	10.1	9.6	0.161	7.1	9.4	9.4	5.3	88	145
OrgSolv82	10.4	9.7	0.26	5.6	8.2	8.5	4.6	88	153
OrgSolv86	10.1	9.5	0.226	6.7	9.1	9.2	5.2	93	155
Acryl10	10.2	9.1	0.164	6.7	10.8	10.6	3.7	89	127
Acryl12	9	9	0.173	5.6	10.1	10	3.1	96	138
Acryl13	9.2	8.9	0.144	6.6	9.7	8.9	5.2	105	92
Acryl1	9.6	9.1	0.184	5.6	10	10	3	88	88
Acryl2	9.6	9	0.185	5.5	10.2	10.1	3	88	132
Acryl22	10	9.2	0.2	5.8	9.5	9.5	3.9	89	105
Acryl24	9.7	9.1	0.132	6.8	9.3	9.4	5.6	98	146
Acryl25	9.2	9.1	0.193	5.7	10	9.9	4.3	91	154
Acryl26	9.2	9.1	0.185	6.8	10.4	10.3	4.8	100	163
Acryl4	9.5	9	0.164	5.7	10	9.9	3.1	89	135
Acryl5	9	9.6	0.131	6.2	10.6	9.6	4.7	98	83
Acryl6	8.9	8.9	0.181	5.8	10.4	10.3	3.6	92	150
Acryl42	9.1	9	0.196	6.1	11	10.8	3	92	109
Acryl8	9.4	9.1	0.19	5.8	9.8	9.7	4	92	150
Acryl9	10	9.2	0.201	6.7	10.5	10.3	4.6	93	148
Acryl14	9.8	9.3	0.157	6.3	10.9	10.7	3.2	88	134
Acryl16	10.3	8.7	0.169	6.9	10.7	10.5	4.6	95	94
Acryl17	9.4	9.1	0.226	6.6	10.3	10.2	3.9	89	146
Acryl30	9.2	9.1	0.173	7.1	10.8	10.6	4.4	96	147
Acryl31	8.8	9.1	0.133	6.4	10.2	10.1	4.2	90	108
Acryl33	9.1	9.1	0.173	6.2	10	9.9	4.6	92	140
Acryl38	9.5	9.1	0.155	6.2	11.2	11	3	88	86
Acryl39	9.7	9	0.163	7.3	10	11.9	3.5	87	87
Acryl18	9	9	0.178	6.3	10	10	4.1	96	141
Acryl19	8.8	9.5	0.171	6.3	10.9	10.7	3.1	92	139
Acryl20	9	9	0.124	6.3	10.5	10.3	3	93	140
aC1	11.5	9.3	0.12	6.8	8.1	8.5	5.4	82	144
aC2	10.6	9.1	0.12	6.7	8.2	8.6	5	92	152
aC3	10.4	9.1	0.1	6.5	8	8.4	5.2	90	150
aC4	10.9	9.1	0.11	6.4	8.1	8.4	. 4.6	90	150
aC5	11.1	9.3	0.13	7.5	8.5	8.8	5.8	88	154
aC7	11	9.2	0.08	7.4	8.3	8.6	5.8	84	154
aC8	11.6	9.2	0.11	7.5	8.5	8.8	5.8	84	148
C13	11.9	9.4	0.1	6.9	8.4	8.7	5.2	80	176
C14	11.4	9.5	0.09	7.4	8.8	9	5	88	166
C16	11.4	9.2	0.1	6.5	8.8	9	3.6	92	176
C18	10.6	9.1	0.1	6.8	8.5	8.8	5.2	90	168

C20	11	9	0.16	6.6	8.4	8.7	4.6	88	168
C15	12.3	9.6	0.1	7.1	9.7	9.6	3.6	88	164
C17	11.3	9.3	0.13	6.8	8.6	8.9	4.6	86	174
C21	11.1	9.3	0.09	6.6	8.6	8.9	4.6	90	192
C22	11.4	9.2	0.13	6.7	8.4	8.6	5	86	176
rC28	10.9	8.3	0.16	7	8.5	8.8	4.9	89	132
rC31	11.1	8.3	0.13	6.3	8.5	8.7	4.5	83	145
rC41	11.3	8.2	0.14	6.6	8.3	8.6	4.7	86	136
rC46	10.7	8.3	0.13	6.6	8	8.4	4.5	83	128
rC47	12.9	8.1	0.16	7.2	8.4	8.6	5.6	74	140
rC48	10.3	8.2	0.13	7.1	8.4	8.7	5.2	93	145
rC51	10.3	8.1	0.11	7	8.5	8.7	4.7	89	135
rC52	10.7	8.3	0.12	7.5	8.1	8.5	5.7	86	149
rC59	10.1	8	0.12	6.8	8.2	8.6	5.1	97	137
rC61	10.7	8.1	0.16	6.9	8.2	8.5	5.1	85	146
rC68	10.6	8.2	0.12	7.1	8.4	8.7	5.6	90	141
rC69	10.4	8.1	0.14	7	8.3	8.6	5.1	92	143
rC71	11	8.2	0.12	7.1	8.1	8.5	5.4	87	128

#### wherein

5

20

- WBC is the white blood cell count (k/ml)
- RBC is the red blood cell count.(M/ml)
- Hb is the hemoglobin content (g/dl)
  - HCT is the hematocrit; volume percentage of erythrocytes (%)
  - MCV is the mean cell volume of erythrocytes  $(fl(\mu^3))$
  - MCH is the mean cell hemoglobin content (pg)
  - MCHC is the mean cell hemoglobin concentration (g/dl)
- PLT is platelet count (k/ml)
  - LYM is the ratio of the number of lymphocytes to the number of leukocytes (%)
  - LYC is the lymphocyte count (k/ml)
  - RDW is red blood cell size distribution width
- RDS is the red blood cell size distribution skewness
  - PCT is the platelet crit; volume percentage of platelets in whole blood (%)
  - MPV is the mean platelet volume  $(fl(\mu^3))$
  - PDW is the platelet size distribution width.
  - PDS is the platelet size distribution skewness.
  - MP is the most frequent platelet volume (fl (μ³))
  - MR is the most frequent red blood cell volume  $(fl(\mu^3))$
  - MW is the most frequent white blood cell volume.  $(fl(\mu^3))$

A test set of the same data was carried out with other patients of the same character as set forth in the preceding tables in order to evaluate the predictive capabilities of the analytic model.

TABLE 2

Test set	WBC	RBC	нв	HCT	MCV	мсн	MCHC	PLT	LYM	LYC
Patient code	WBC	RBC	nB	nc1	MCV	MCH	Hene	1 11	J	
A2C3	10.4	5.48	15.8	45.9	83.9	28.8	34.4	242	20	1
A3C3	8.2	4.93	14.8	43.2	87.7	30	34.3	253	30.2	2.5
A3K3	6.9	5.33	16.1	44	82.9	30.3	36.6	139	37.4	2.6
A213	5.8	5.82	18	50	85.9	31	36.1	158	27.8	1.6
Amlgm5	4.6	4.43	15.7	42.2	95.3	35.5	37.3	166	33.8	1.6
Amigm12	7.8	4.36	13.9	38	87.1	31.9	36.6	136	30.1	2.4
Amlgm21	3.9	5.3	17.1	48.2	90.8	32.4	35.6	127	42.5	1.6
Amlgm22	4.7	4.83	17.2	48.2	99.7	35.5	35.5	173	40.9	1.9
OrgSolv75	6.5	4.55	15.7	41.1	90.3	34.5	38.2	286	33.9	2.2
OrgSolv97	3.2	5.3	16.5	48.9	91.4	30.9	33.8	329	55	1.8
OrgSolv69	5.9	4.85	16.6	41.5	85.5	34.1	39.9	273	42.3	2.6
OrgSolv83	3.8	4.8	15.6	46.8	96.8	32.3	33.4	270	41.3	1.6
OrgSolv87	6.1	5.6	18.6	55.1	97.6	32.8	33.6	287	31.2	1.9
Acryl3	4.1	4.45	13.5	41	92	30.4	33.1	30.7	41.9	1.7
Acryl21	4.4	4.59	14.3	43.3	94.4	31.5	33.1	307	43.2	2
Acryl43	5.1	4.38	13.4	40.5	92.5	30.6	33.1	367	26.7	1.4
Acryl7	5.4	5.13	15	46.4	90.4	29.4	32.5	310	42.8	2.3
АсгуИ0	4.3	4.17	12.7	39.3	94	30.5	32.4	210	39.7	1.7
aC6	5.8	4.78	12.4	40.4	84.6	25.9	30.7	199	35.7	2.1
C19	7.8	4.61	17.8	43.2	93.8	38.7	41.2	177	25.6	
C12	5.1	5.04	17.3	43.4	86.3	34.3	39.8	150	41.4	2.1
rC73	5.2	5.28	15.5	42.9	81.3	29.3	36.1	181	36.3	1.9
rC81	5.1	5.18	15.5	43.3	83.5	30	35.9	222	36.4	1.8
Patient code	RDW	RDS	PCT	MPV	PDW	PDS	MP	MR	MW	
A2C3	10.6	8.1	0.16	6.5	8.2	8.5	5	93	134	
A3C3	11.2	8.1	0.171	6.8	8.2	8.5	5.4	98	151	
A3K3	8.7	8	0.096	6.9	8.7	8.9	4.4	93	147	
A2I3	10.9	8.2	0.11	7	8.4	8.7	5.6	95	139	
Amlgm5	10	9.1	0.11	6.8	8.5	8.8	5.8	94	168	
Amlgm12	11.9	9.3	0.1	7.2	8.4	8.7	6.2	84	15B	
Amigm21	. 10	9.2	0.09	6.8	8.3	8.7	5.4	88	180	
Amlgm22	10.8	9.1	0.12	7.2	8.5	8.8	6.2	100	182	
OrgSolv75	9.2	8.6	0.166	5.8	8.7	8.9	4.3	94 84	138 94	
OrgSolv97	11.7	9.8	0.248	7.5	9.5	9.6	6.2	85	85	
OrgSolv69	10.7	9	0.174	6.4	9.5	9.6			105	
OrgSolv83	10.5	9.5	0.177	6.5	9.2	9.2	5.3	91 92	105	
OrgSotv87	10.3	9.5	0.179	6.2	9.6	9.6	4.9	92	84	
Acryl3	9.4	8.9	0.174	5.6	10.6	10.5	3.9	93	83	
Acryl21	10	9.2	0.197	6.4	10.4	10.3		93	151	
Acryl43	8.7	9	0.222	6	10.1	10	3.1	89	89	
Acryl7	10.1	9.1	0.195	6.2	10.5	10.3	3.1	94	85	
Acryl40	8.7	9	0.13	6.2	10.6	10.4	5.6	82	150	
aC6	10.8	9.2	0.15	7.3	8.4	8.7 8.8	5.6	92	166	
C19	11.2	9.2	0.12	6.8	8.6	8.8	5.2	82	162	
C12	12	9.5	0.11	7.5	8.8	8.9	4.7	90		
rC73	11	8.1	0.11	6.1	8.1			90	135	
rC81	10.4	8	0.14	6.4	8.3	8.6	4.4	92	135	

The above blood cell analyses are copied from Excel® into Matlab® for Windows™, version 4.0, sold by MathWorks, Inc., Natick, Massachusetts, and autoscaled there using the following routine in Matlab® for Windows™. The table of data is designated "x" and the code for autoscaling is —

Function i=autosc(x)
[m,n]=size(x);
i=x./(ones(m,1)\*std(x));
i=i-ones(m,1)\*mean(i);
end

5

10

15

The next step is the orthogonalization of the matrix. The autoscaled data is subjected to the singular value decomposition,

[U,S,V]=svd(autoscaled data matrix)

which is a command in Matlab® for Window.™ In order to produce the input data, a multiplication U times (\*) S is carried out in Matlab® for Window.™ The orthogonalized matrix comprises the input data ("q") for the artificial neural network.

The following table 3 is the output data corresponding to the input data ("q").

TABLE 3

Patient code	Heart	Amalgam	OrgSolv	Acryl	Controls
A1A3	1	0	0	0	0
A1B3	1	0	0	0	0
A1G3	1	0	0	0	0
A2A3	1	0	0	0	0
A2D3	1	0	0	0	0
A2E3	1	0	0	0	0
A2G3	1	0	0	0	0
A2H3	1	0	0	0	0
A2J3	1	0	0	0	0
A2K3	1	0	0	0	0
A2L3	1	0	0	0	0
A3B3	1	0	0	0	0
A3D3	1	0	0	0	0
A3E3	1	0	0	0	0
A3F3	1	0	0	0 •	0
A3G3	1	0	0	0	0
A3H3	1	0	0	0	0
A3I3	1	0	0	0	0
A3M3	1	0	0	0	0
Amlgm1	0	1	0	0	0
Amlgm2	0	1	0	0	0
Amlgm3	0	1	0	0	0
Amlgm4	0	1	0	0	0
Amlgm6	0	1	0	0	0
Amlgm7	0	1	0	0	0
Amlgm8	0	1	0	0	0
Amigm9	0	1	0	0	0
Amlgm10	0	1	0	0	0
Amlgm11	0	1	0	0	0 .
Amlgm13	0	1	0	0	0
Amlgm14	0	1	0	0	0
Amlgm15	0	1	0	0	0
Amlgm16	0	1	0	0	0
Amlgm17	0	1	0	0	0

		4	0	0	0
Amlgm20	0	1	0	0	0
Amlgm23	0		0	0	0
Amlgm24	0	1	0	0	0
Amlgm25	0	1		0	0
OrgSolv73	0	0	1		0
OrgSolv74	0	0	11	0	
OrgSolv98	0	0	1	0	0
OrgSolv99	0	0	1	0	0
OrgSolv67	0	0	1	0	0
OrgSolv68	0	0	1	0	0
OrgSolv84	0	0	1	0	0
OrgSolv85	0	0	1	0	0
OrgSolv94	0	0	1	0	0
OrgSolv95	0	0	1	0	0
OrgSolv96	0	0	1	0	0
OrgSolv70	0	0	1	0	0
OrgSolv71	0	0	1	0	0 .
OrgSolv72	0	0	1	0	0
OrgSolv64	0	0	1	0	0
OrgSolv65	0	0	1	0	0
OrgSolv66	0	0	1	0	0
OrgSolv91	0	0	1	0	0
OrgSolv92	0	0	1	0	0
	0	Ö	1	0	0
OrgSolv93	0	0	<del>- i</del>	Ö	0
OrgSolv88		0	1	0	ō
OrgSolv89	0		1	0	0
OrgSolv90		0	1	0	0
OrgSolv61	0	0		0	0
OrgSolv62	0	0	1	0	0
OrgSolv63	0	0	1		0
OrgSolv76	0	0	11	0	
OrgSolv77	0	0	11	0	0
OrgSolv78	0	0	11	0	
OrgSolv60	0	0	11	0	0
OrgSolv79	0	0	1 .	0	0
OrgSolv80	0	0	1	0	0
OrgSolv81	0	0	1	0	0
OrgSolv82	0	0	1	0	0
OrgSolv86	00	0	1	0	0
Acryl10	0	0	0	1	0
Acryl12	0	0	0	1	0
Acryl13	0	0	0	1	0
Acryl1	0	0	0	1	0
Acryl2	0	0	0	1	0
Acryl22	0	0	0	1	0
Acryl24	0	0	0	1	0
Acryl25	0	0	0	1	0
Acryl26	0	0	0	1	0
Acryl4	0	0	0	1	0
Acryl5	0	0	0	1	0
Acryl6	0	0	0	1	0
Acryl42	0	0	0	1	0
Acryl8_	Ō	0	0	1	0
Acryl9	0	0	0	1	0
Acryl14	0	0	ō	1	0
Acryl16	0	0	ŏ	1	0
Acryl17	0	0	0	1	0
	0	0	0	1	0
Acryl30	0	0	0	1	0
Acryl31	<u> </u>			<u> </u>	<u> </u>

Acryl38	Acryl33	0	0	0	1	0
Acryl39						
Acryl18						
Acryl19         0         0         0         1         0           Acryl20         0         0         0         1         0           aC1         0         0         0         0         1           aC2         0         0         0         0         1           aC3         0         0         0         0         1           aC3         0         0         0         0         1           aC4         0         0         0         0         1           aC5         0         0         0         0         1           aC7         0         0         0         0         1           aC8         0         0         0         0         1           aC8         0         0         0         0         1           aC13         0         0         0         0         1           C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1						
Acryl20						
aC1         0         0         0         0         1           aC2         0         0         0         0         1           aC3         0         0         0         0         1           aC4         0         0         0         0         1           aC5         0         0         0         0         1           aC5         0         0         0         0         1           aC7         0         0         0         0         1           aC8         0         0         0         0         1           aC8         0         0         0         0         1           C13         0         0         0         0         1           C14         0         0         0         0         1           C14         0         0         0         0         1           C14         0         0         0         0         1           C18         0         0         0         0         1           C19         0         0         0         0         1 <t< td=""><td></td><td>1</td><td></td><td></td><td></td><td></td></t<>		1				
aC2         0         0         0         0         1           aC3         0         0         0         0         1           aC4         0         0         0         0         1           aC5         0         0         0         0         1           aC7         0         0         0         0         1           aC8         0         0         0         0         1           C13         0         0         0         0         1           C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1 <t< td=""><td></td><td>1</td><td></td><td></td><td></td><td></td></t<>		1				
aC3         0         0         0         0         1           aC4         0         0         0         0         0         1           aC5         0         0         0         0         0         1           aC7         0         0         0         0         0         1           aC8         0         0         0         0         0         1           C13         0         0         0         0         0         1           C14         0         0         0         0         0         1           C14         0         0         0         0         0         1           C16         0         0         0         0         0         1           C18         0         0         0         0         0         1           C20         0         0         0         0         1         1           C15         0         0         0         0         1         1           C17         0         0         0         0         0         1         1           C21		1				
aC4         0         0         0         0         1           aC5         0         0         0         0         0         1           aC7         0         0         0         0         0         1           aC8         0         0         0         0         0         1           C13         0         0         0         0         0         1           C14         0         0         0         0         0         1           C16         0         0         0         0         0         1           C18         0         0         0         0         0         1           C20         0         0         0         0         0         1           C20         0         0         0         0         1         1           C15         0         0         0         0         1         1           C21         0         0         0         0         1         1           C21         0         0         0         0         0         1         1           C22						
aC5         0         0         0         0         1           aC7         0         0         0         0         1           aC8         0         0         0         0         0         1           C13         0         0         0         0         0         1           C14         0         0         0         0         0         1           C16         0         0         0         0         0         1           C18         0         0         0         0         0         1           C18         0         0         0         0         0         1           C20         0         0         0         0         0         1           C20         0         0         0         0         1         1           C15         0         0         0         0         0         1         1           C21         0         0         0         0         0         1         1         1         1         1         1         1         1         1         1         1         1         1 <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>						
aC7         0         0         0         0         1           aC8         0         0         0         0         1           C13         0         0         0         0         1           C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1						
aC8         0         0         0         0         1           C13         0         0         0         0         1           C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           C22         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1						
C13         0         0         0         0         1           C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           C23         0         0         0         0         1           C24         0         0         0         0         1           C31         0         0         0         0         1           C44         0         0         0         0         1 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td></t<>						
C14         0         0         0         0         1           C16         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C15         0         0         0         0         0           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           C28         0         0         0         0         1           C28         0         0         0         0         1           C28         0         0         0         0         1           C31         0         0         0         0         1           C41         0         0         0         0         1 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td></t<>						
C16         0         0         0         0         1           C18         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           C28         0         0         0         0         1           C28         0         0         0         0         1           C28         0         0         0         0         1           C31         0         0         0         0         1           C41         0         0         0         0         1           C44         0         0         0         0         1           C46         0         0         0         0         1           rC47         0         0         0         0         1           <						
C18         0         0         0         0         1           C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           C28         0         0         0         0         1           C31         0         0         0         0         1           C41         0         0         0         0         1           C41         0         0         0         0         1           C46         0         0         0         0         1           C47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1		0	0			1 .
C20         0         0         0         0         1           C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           rC28         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC52         0         0         0         0         1           rC61         0         0         0         0         1	C16	I	0		0	1
C15         0         0         0         0         1           C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           rC28         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC69         0         0         0         0         1	C18	0	0	0	0	1
C17         0         0         0         0         1           C21         0         0         0         0         1           C22         0         0         0         0         1           rC28         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC69         0         0         0         0         1	C20	0	0	0	0	1
C21         0         0         0         0         1           C22         0         0         0         0         1           rC28         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1	C15	0	0	0	0	1
C22         0         0         0         0         1           rC28         0         0         0         0         1           rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1	C17	0	0	0	0	1
rC28         0         0         0         0         1           rC31         0         0         0         0         0         1           rC41         0         0         0         0         0         1           rC46         0         0         0         0         1         1           rC47         0         0         0         0         1	C21	0	0	0	0	1
rC31         0         0         0         0         1           rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1	C22	0	0	0	0	1
rC41         0         0         0         0         1           rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         0	rC28	0	0	0	0	1
fC41         0         0         0         0         1           fC46         0         0         0         0         1           fC47         0         0         0         0         1           fC48         0         0         0         0         1           fC51         0         0         0         0         1           fC52         0         0         0         0         1           fC59         0         0         0         0         1           fC61         0         0         0         0         1           fC68         0         0         0         0         1           fC69         0         0         0         0         0	rC31	0	0	0	0	1
rC46         0         0         0         0         1           rC47         0         0         0         0         1           rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1	rC41	0	0	0	0	1
fC47         0         0         0         0         1           fC48         0         0         0         0         1           fC51         0         0         0         0         1           fC52         0         0         0         0         1           fC59         0         0         0         0         1           fC61         0         0         0         0         1           fC68         0         0         0         0         1           fC69         0         0         0         0         1		0	0	0	0	1
rC48         0         0         0         0         1           rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1		<u> </u>	0	0	0	1
rC51         0         0         0         0         1           rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1				0	0	1
rC52         0         0         0         0         1           rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1		<del></del>		0	0	1
rC59         0         0         0         0         1           rC61         0         0         0         0         1           rC68         0         0         0         0         1           rC69         0         0         0         0         1				0	0	1
rC61 0 0 0 0 1 1 rC68 0 0 0 0 0 1 rC69 0 0 0 0 1				0	0	1
rC68 0 0 0 0 0 1 1 rC69 0 0 0 0 1						1
rC69 0 0 0 0 1		<u> </u>				1
1003						1

		TABLE	4					
Test Set Data:								
	1	0	0	0	0			
	1	0_	0	0	0			
	1	0	0	0	0			
	1	0 ·	0	0	0			
	0	1	0	0	0			
	0	1	0	0	0			
	0	1	0	0	0			
	0	1	0	0	0			
	0	0	1	0	0			
	0	0	1	0	0			
	0	0	1	0	0			
	0	0	1	0	0			
	0	0	1	0	0			
	0	0	0	1	0			
	0	0	0	1	0			
	0	0	0	1	0			
	0	0	0	1	0			
	0	0	0	1	0			
	0	0	0	0	1			
	0	0	0	0	1			

0	0	0	0	1
 0	0	0	0	1
 Ó	0	0	0	1

#### The Artificial Neural Network

The artificial neural network is configured using Matlab® for Windows™, according to the number of input variables, which in the above case (Table 1) is 19, providing the number of input neurons.

The following program<sup>2</sup> is generated on a computer by Matlab® for Windows<sup>74</sup> to configure the artificial neural network by inputting the first line with function and initff data:

$$[w1,b1,w2,b2] = initff(p,s1,f1,T,f2)$$

The arguments in the above equation for the above illustration are as follows:

- w1 is the weight matrix for the hidden layer;
- b1 is the vector of biases for the hidden layer;
- w2 is the weight matrix for the output layer;
- b2 is the vector of biases for the output layer;
- p is the matrix of minimum and maximum values for input matrix;
- sl is the number of neurons in the hidden layer which is 40 in this case;
  - f1 is the transfer function for the hidden layer and in this case is called the hyperbolic tangent sigmoid function:
- T is matrix of target vectors as set forth in Table 3 above; and
  - f2 is the transfer function for the output layer, which in this case it is the logistic sigmoid function.

%INITFF Initialize feed-forward network up to 3 layers.

30 8

10

15

20

25

 $\{W1,B1,\ldots\} =$ 

INITFF(P,S1,'F1',...,Sn,'Fn')

Revision: 1.1 Date: 1994/01/11 16:24:46.

 $<sup>^{2}</sup>$  Written by Mark Beale, 12-15-93 and Copyright © 1992-94 by the MathWorks,

WO 97/17891

PCT/US96/18254

20

```
P - Rx2 matrix of input vectors.
                                Si - Size of ith layer.
       ģ
                                Fi - Transfer function of the ith
       layer (string).
5
                                Returns:
                                         Wi - Weight matrix of the
       ith layer.
                                         Bi - Bias (column) vector of
       the ith layer.
                                 [W1,B1,\ldots] =
10
       INITFF(P,S1,'F1',...,T,'Fn')
                                 T - SnxQ matrix of target vectors.
                                  Returns weights and biases.
                                  IMPORTANT: Each ith row of P must
15
       contain expected
                                   min and max values for the ith in-
       put.
       g,
                                   EXAMPLE: [W1,b1,W2,b2] = initff([0
20
       10; -5 5],5,'tansig',3,'pure
                                               p = \{4; -1\};
                                    a =
       simuff(p,W1,b1,'tansig',W2,b2,'purelin')
                                    See also NNINIT, BACKPROP,
25
       SIMUFF, TRAINBPX, TRAINLM
                                       if all([3 5 7] ~= nar-
       gin), error('Wrong number of input arguments'), end
        [R,Q] = size(p);
        if (Q < 2), error('P must contain at least two elements in
30
        each row.'), end
        if nargin == 3
          [S,Q] = size(s1);
          if max(S,Q) > 1, s1 = S; end
          [w1,b1] = feval(feval(f1,'init'),s1,p);
35
        elseif nargin == 5
```

```
21
```

[w1,b1] = feval(feval(f1,'init'),s1,p); x = ones(s1,1) \* feval(f1,'output');[S,Q] = size(s2);if max(S,Q) > 1, s2 = S; end [w2,b2] = feval(feval(f2,'init'),s2,x);5 elseif nargin == 7 [w1,b1] = feval(feval(f1,'init'),s1,p); x = ones(s1,1) \* feval(f1,'output');[w2,b2] = feval(feval(f2,'init'),s2,x);10 x = ones(s2,1) \* feval(f2,'output');[S,Q] = size(s3);if max(S,Q) > 1, s3 = S; end [w3,b3] = feval(feval(f3,'init'),s3,x);15 end

> The next step involves training the artificial neural network. Again, the computation is carried out in a computer using Matlab® for Windows™ through use of the program called "TrainBPX.M". 3 In this case, input is made to the equation

20

25

30

[W1,B1,W2,B2,...,TE,TR] =TRAINBPX(W1, B1, F1, W2, B2, F2, ..., P, T, TP)

wherein

- W1 on the left side of the equation is the initial weight matrix for the hidden layer and on the right side of the equation, the optimal weight matrix;
- B1 on the left side of the equation is the initial vector of biases for the hidden layer and on the right of the equation is the optimal vector of biases for the hidden layer;
- W2 on the left side is the initial weight matrix for the output layer and on the right side is the optimal weight matrix for the output layer;

The program is indicated to be copyrighted by Mark Beale and by the MathWorks, Inc.

PCT/US96/18254 WO 97/17891 22

• B2 on the left side is the vector of biases for the output layer and on the right side is the optimum vector of biases for the output layer;

- P is the matrix of input vectors as recited in Table 1;
- TE is the actual number of training passes or epics (in this case 5000);
  - TR is the training record which is the stored error of pass through the network;
  - F1 is the transfer function for the hidden layer and in this case is called the hyperbolic tangent sigmoid
    - T is matrix of target vectors as set forth in Table 3 above;
    - F2 is the transfer function for the output layer, which in this case it is the logistic sigmoid function; and
    - TP is an optional vector of training parameters. The Matlab® for Windows™ TrainBPX.M screen appears as:

"function [a,b,c,d,e,f,g,h]

20 trainbpx(i,j,k,l,m,n,o,p,q,r,s,t) %TRAINBPX Train feed-forward network with fast backpropagation.

5

10

15

- TRAINBPX can be called with 1, 2, or 3 sets of weights
- and biases to train up to 3 layer feed-forward networks. 25

[W1, B1, W2, B2, ..., TE, TR]

TRAINBPX (W1, B1, F1, W2, B2, F2, ..., P, T, TP)

- Wi Weight matrix for the ith layer.
- Bi Bias vector for the ith layer. 30
  - Fi Transfer function (string) for the ith layer.
  - P RxQ matrix of input vectors.
  - T SxQ matrix of target vectors.
  - TP Training parameters (optional).
- Returns new weights and biases and 35 å
  - Wi new weights.
  - Bi new biases.

23

```
TE - the actual number of epochs trained.
               TR - training record: [row of errors]
            Training parameters are:
               TP(1) - Epochs between updating display, default = 25.
5
               TP(2) - Maximum number of epochs to train, default =
      1000.
              TP(3) - Sum-squared error goal, default = 0.02.
              TP(4) - Learning rate, 0.01.
              TP(5) - Learning rate increase, default = 1.05.
10
              TP(6) - Learning rate decrease, default = 0.7.
               TP(7) - Momentum constant, default = 0.9.
               TP(8) - Maximum error ratio, default = 1.04.
            Missing parameters and NaN's are replaced with defaults.
15
             See also NNTRAIN, BACKPROP, INITFF, SIMFF, TRAINBP,
       TRAINLM.
          % Mark Beale, 9-2-92
          % Revised 12-15-93, MB.
          % Copyright (c) 1992-94 by the MathWorks, Inc.
20
          % $Revision: 1.1 $ $Date: 1994/01/11 16:30:06 $
       if all([5 6 8 9 11 12] ~= nargin),error('Wrong number of in-
       put arguments'), end
       if nargin == 5
25
         [a,b,c,d] = tbpxl(i,j,k,l,m);
       elseif nargin == 6
         [a,b,c,d] = tbpxl(i,j,k,l,m,n);
       elseif nargin == 8
         [a,b,c,d,e,f] = tbpx2(i,j,k,1,m,n,o,p);
       elseif nargin == 9
30
         [a,b,c,d,e,f] = tbpx2(i,j,k,l,m,n,o,p,q);
       elseif nargin == 11
         [a,b,c,d,e,f,g,h] = tbpx3(i,j,k,l,m,n,o,p,q,r,s);
       elseif nargin == 12
         [a,b,c,d,e,f,g,h] = tbpx3(i,j,k,l,m,n,o,p,q,r,s,t);
35
       end"
```

The next step involves using the trained artificial neural network and the analytic model thereof to obtain a diagnosis for a patient or patients that have not previously been diagnosed. This is done by the following sequence:

- withdrawing from a new patient that has not been diagnosed for such disorder, a sample or samples of a kind taken from said reference patients and subjecting said diagnostic sample to said clinical testing to obtain new numerical data;
- automatically scaling in the computer the new data to the data derived from the reference patients;
- feeding the scaled new data to the trained artificial neural network and the analytic model thereof, and
- automatically obtaining a diagnosis as to the new patient with respect to such disorders.

Again, the computation is carried out on a computer using Matlab® for Windows™ through use of the program called "Simuff..M". In this case, input is made to the equation —

[A] = SIMUFF(P,W1,B1,'F1',W2,B2,'F2')

#### 20 wherein

5

10

15

25

30

- A is the result from the layer of output neurons, i.e. the prediction;
- W1 is the optimal weight matrix;
- B1 is the optimal vector of biases for the hidden layer;
- W2 is the optimal weight matrix for the output layer;
- B2 is the optimum vector of biases for the output layer;
- P is the matrix of input vectors obtained from undiagnosed patients;
- F1 is the transfer function for the hidden layer and in this case is called the hyperbolic tangent sigmoid function;

The program is indicated to be copyrighted by Mark Beale and by the MathWorks, Inc.

PCT/US96/18254 WO 97/17891 25

```
    F2 is the transfer function for the output layer, which

              in this case it is the logistic sigmoid function
             Simuff..M appears on the computer screen as:
       "function [a1,a2,a3] = simuff(p,w1,b1,f1,w2,b2,f2,w3,b3,f3)
 5
       %SIMUFF Simulate feed-forward network.
             SIMUFF will simulate networks with up to 3 layers.
        ક
             SIMUFF(P,W1,B1,'F1',...,Wn,Bn,'Fn')
10
        ş
               P - Matrix of input (column) vectors.
               Wi - Weight matrix of the ith layer.
               Bi - Bias (column) vector of the ith layer.
        욯
               Fi - Transfer function of the ith layer (string).
             Returns output of nth layer.
15
             [A1, A2, ...] = SIMUFF(P, W1, B1, 'F1', ..., Wn, Bn, 'Fn')
             Returns:
               Ai - Output of the ith layer.
                                               initff([0
                                                              10;
                                                                     -5
20
             EXAMPLE:
                         [w1,b1,w2,b2] =
       5], 3, 'tansig', 2, 'purelin');
                      p = [2; -3];
                      a = simuff(p,w1,b1,'tansig',w2,b2,'purelin')
             See also NNSIM, BACKPROP, INITFF, TRAINBPX, TRAINLM.
25
       if all([4 7 10] ~= nargin),error('Wrong number of input argu-
       ments'), end
       if nargin == 4
30
         al = feval(f1, wl*p, b1);
       elseif nargin == 7
         al = feval(fl, wl*p, bl);
         a2 = feval(f2, w2*a1, b2);
         if nargout <= 1, a1 = a2; end
35
       elseif nargin == 10
         al = feval(fl, wl*p, bl);
         a2 = feval(f2, w2*a1, b2);
         a3 = feval(f3, w3*a2, b3);
```

26

if nargout <= 1, a1 = a3; end
end</pre>

5

10

15

20

25

30

35

The validity of the above described artificial neural network analytic model can be tested and quantified using a standard procedure in multivariate modeling, namely that of an independent test set comprising data from a group of patients that have been previously diagnosed, described in table 4, and have had samples withdrawn of a kind taken from said reference patients. The diagnostic sample is subjected to said clinical testing to

- obtain new numerical data, described in table 2,
- automatically scaling in the computer the new data to the data derived from the reference patients,
- feeding the scaled new data to the trained artificial neural network and the analytic model thereof, and
- automatically obtaining a diagnosis as to the test-set patients with respect to such disorders.

The diagnosis predicted by the artificial neural network can be compared to the diagnosis previously obtained and one of a number of different measures of the validity of the artificial neural network and the analytic model thereof can be calculated as follows: Classification Success (%);

CS = 100  $\star$  (wrong classifications/classifications) A typical value of CS for the artificial neural network is 92%.

As a final step in the process an approximate probability level, associated to the diagnosis of disorder via the above described artificial neural network, can be calculated electronically in the computer using a variant of the Students t-test to obtain the probability level.

The calculation of the probability level is done in a three step process. The first step involves building a PLS-Discriminant Analysis model (M. Sjöström, S. Wold and B. Söderström, PLS Discriminant Plots, in "Pattern Recognition in Practice II" (1986), E. S. Gelsema and L. N. Kanal, Eds., Elsevier Science Publishers (North-Holland), pages 461-470) for each disorder group, using the aforementioned reference groups. The second step uses the above described artificial

neural network analytic model to select the matching disorder for the unknown and undiagnosed patient and the associated PLS-Discriminant Analysis model. The third and final step uses the selected PLS-Discriminant Analysis model and the above mentioned Students t-test to ascertain the probability level of the diagnosis.

In the first step, a PLS-Discriminant Analysis model for each reference group is developed, separately for each reference group. Again, the computation is carried out using Matlab® for Windows<sup>™</sup> through use of the program called "PLS.M". In this case, input is made to the equation;

[P,Q,W,T,U,b,ssqdif] = pls(x,y,lv)

#### wherein

WO 97/17891

5

10

- P is the matrix of x-loadings;
- Q is the matrix of y-loadings;
  - W is the matrix of weights;
  - T is the matrix of x-scores;
  - U is the matrix of y-scores;
  - b is the vector of regression coefficients between x and y;
- ssqdif is the variance contained in each latent variable;
  - x is the matrix of input vectors, comprised of data from one reference group and one group of healthy patients, all as recited in table 5;
  - y is the matrix of target vectors as in table 6;
- lv is the total number of latent variables to be calculated function [P,Q,W,T,U,b,ssqdif] = pls(x,y,lv)
  - % This PLS program calculates the P, Q, W, T, and U matrices
  - for the given x and y blocks and the number of latent
  - % variables specified. It also outputs the fraction of
- 30 % variance used in the x and y matrices.
  - % Input format is as follows:
  - % function [P,Q,W,T,U,b,ssqdif] = pls(x,y,lv);

[mx,nx] = size(x);

 $<sup>^{5}</sup>$  The computer screen indicates that the program is the copyright of Barry M. Wise, 1991.

WO 97/17891 28 PCT/US96/18254

```
[my, ny] = size(y);
       if nx < lv
       error('No. of LVs must be <= no. of x-block variables')
 5
       P = zeros(nx, lv);
       Q = zeros(ny, lv);
       W = zeros(nx, lv);
       T = zeros(mx, lv);
       U = zeros(my, lv);
       b = zeros(1, lv);
10
       ssq = zeros(lv, 2);
       ssqx = 0;
       for I = 1:nx
         ssqx = ssqx + sum(x(:,I).^2);
15
       end
       ssqy = sum(sum(y.^2)');
       for I = 1:1v
          [p,q,w,t,u] = plsnipal(x,y);
         b(1,I) = u'*t/(t'*t);
         x = x - t*p';
20
         y = y - b(1, I) *t*q';
         ssq(I,1) = (sum(sum(x.^2)'))*100/ssqx;
          ssq(I,2) = (sum(sum(y.^2)'))*100/ssqy;
         T(:,I) = t(:,1);
         U(:,I) = u(:,1);
25
          P(:,I) = p(:,1);
         W(:,I) = W(:,1);
          Q(:,I) = q(:,1);
        end
        ssqdif = zeros(lv,2);
30
        ssqdif(1,1) = 100 - ssq(1,1);
        ssqdif(1,2) = 100 - ssq(1,2);
        for I = 2:lv
          for j = 1:2
            ssqdif(I,j) = -ssq(I,j) + ssq(I-1,j);
35
          end
        end
        disp(' ')
```

29

```
disp(' Percent Variance Captured by PLS Model')
disp(' ')
disp(' ----X-Block------')
disp(' LV# This LV Total This LV Total ')
disp([(1:lv)' ssqdif(:,1) cumsum(ssqdif(:,1)) ssqdif(:,2)
cumsum(ssqdif(:,2))])
```

Table 5

Reference	w	R	нв	Н	М	M	М	Р	L	L	R	R	Р	М	P	Р	MP	MR	MW
Patients	В	В		С	С	С	С	L	Y	Y	D	D	С	P	D	D			( )
	С	С		T	V	н	н	Т	М	С	W	s	Т	l v	w	S			į
							С												l
Amigm1	12.5	4.35	17.4	45.1	103.7	40.1	38.7	193	22.5	2.8	11.1	9.2	0.13	6.6	8.6	8.8	5.3	102	176
Amigm2	5.9	4.2	13.4	39.7	94.4	31.9	33.8	135	38	2.3	10.8	9.3	0.1	7.7	9	9.2	5.8	92	154
Amlgm3	8	4.24	13.7	40.6	95.6	32.2	33.7	150	33	2.6	10.8	9.1	0.12	7.7	8.7	8.9	6.5	96	158
Amlgm4	6.7	5.07	19.1	46.5	96.6	37.7	41.2	191	37.1	2.5	10.7	9.2	0.12	6.5	8.6	8.8	5	90	174
Amigm6	6	4.25	14.4	39.6	93.2	34	36.4	112	35.1	2.1	9.8	9.2	0.08	7.4	9	9.2	6.2	92	170
Amigm7	7	4.55	15.5	41.9	92.1	34.1	37	149	26.8	1.9	9.8	9.1	0.11	7.6	8.7	8.8	6.2	92	176
Amigm8	14.2	4.81	17	45.6	94.9	35.4	37.4	169	21.5	3.1	10.5	9.4	0.12	7	8.6	8.8	5.8	92	172
Amlgm9	8.3	4.35	15.3	41.7	95.7	35.1	36.6	159	38.6	3.2	10.6	9.4	0.11	7.2	8.4	8.6	5.6	92	166
Amigm10	8.2	5.44	18	49.3	90.7	33.1	36.5	147	33.3	2.7	11.1	9.3	0.11	7.2	8.6	8.8	6.4	88	168
Amigm11	10.6	4.58	16.3	44	96.1	35.7	37.2	166	33.3	3.5	10.1	9.2	0.12	7.1	8.5	8.8	6.2	94	162
Amlgm13	6.2	4.38	14.3	39.1	89.3	32.6	36.5	169	41.1	2.5	11.8	9.4	0.13	7.4	8.6	8.8	6.2	86	160
Amlgm14	6.1	4.86	17.1	46.2	95.3	35.2	36.9	146	40.9	2.5	10.3	9.2	0.11	7.2	8.8	9	5.2	92	166
Amigm15	8.8	4.85	17.3	48.8	100.5	35.7	35.5	170	28.4	2.5	10.7	9.2	0.11	6.7	8.5	8.8	5.2	98	166
Amlgm16	5.8	3.87	13.3	37.7	97.2	34.4	35.4	158	36.4	2.1	9.8	9.1	0.11	7.1	8.9	9.1	6.4	96	162
Amlgm17	5.6	5.21	16.5	48.8	93.8	31.8	33.9	221	25.6	1.4	11.2	9.3	0.15	6.8	8.9	9.1	4.6	90	160
Amlgm20	5.6	4.89	15.7	44	89.9	32.1	35.7	163	32.4	1.9	10.6	9.2	0.12	7.2	8.4	8.7	6.4	88	160
Amigm23	5.1	4.15	11.7	38.4	92.6	28.1	30.4	201	35.5	1.8	10.8	9.2	0.14	7.2	8.4	8.7	5.4	90	152
Amlgm24	7.2	5.13	13.2	43.8	85.5	25.6	30	240	20.1	1.4	11.9	9.3	0.17	7.1	8.5	8.7	5.8	84	154
Amigm25	6.8	4.41	12	38.6	87.6	27.2	31	233	28.6	1.9	10.8	9.2	0.17	7.2	7.9	8.3	6	86	150
Healthy F	atien	ts																	
aC1	5.5	4.53	11.5	37.8	83.6	25.4	30.4	180	35.7	2	11.5	9.3	0.12	6.8	8.1	8.5	5.4	82	144
aC2	5.9	4.38	12.5	40.7	92.9	28.5	30.7	175	29.9	1.7	10.6	9.1	0.12	6.7	8.2	8.6	5	92	152
aC3	3.9	4.59	12.8	41.4	90.2	27.9	30.9	153	41.7	1.6	10.4	9.1	0.1	6.5	8	8.4	5.2	90	150
aC4	6.1	5.04	14.4	46	91.3	28.6	31.3	175	30.7	1.9	10.9	9.1	0.11	6.4	8.1	8.4	4.6	90	150
aC5	6.3	4.94	14.1	44.9	90.8	28.5	31.4	178	39.7	2.5	11.1	9.3	0.13	7.5	8.5	8.8	5.8	88	154
aC7	2.9	4.74	12.3	40.7	85.6	26	30.3	108	45.9	1.3	11	9.2	0.08	7.4	8.3	8.6	5.8	84	154
aC8	5.2	5.27	14.4	45.7	86.8	27.3	31.9	143	34	1.8	11.6	9.2	0.11	7.5	8.5	8.8	5.8	84	148
C13	4.3	5.09	16.7	43.2	84.9	32.8	38.7	147	34	1.4	11.9	9.4	0.1	6.9	8.4	8.7	5.2	80	176
C14	5.4	4.73	17.9	43.9	93	37.9	40.8	127	36.9	2	11.4	9.5	0.09	7.4	8.8	9	5	88	166
C16	6.1	4.78	17.2	43.4	94.2	35.9	38.2	153	38.1	2.3	11.4	9.2	0.1	6.5	8.8	9	3.6	92	176
C18	7.4	4.88	18.6	45	92.2	38.1	41.3	145	30.B	2.3	10.6	9.1	0.1	6.8	8.5	8.8	5.2	90	168
C20	7	5.25	19.2	47	89.4	36.5	40.8	241	26.9	1.9	11	9	0.16	6.6	8.4	8.7	4.6	88	168
C15	4.8	4.86	18.3	45.6	93.8	37.8	40.2	144	28.3	1.4	12.3	9.6	0.1	7.1	9.7	9.6	3.6	88	164
C17	6.7	5.44	19.6	47.8	87.9	36.1	41.1	193	37.4	2.5	11.3	9.3	0.13	6.8	8.6	8.9	4.6	86	174
C21	6.9	4.82	18.2	44.6	92.5	37.8	40.9	138	23	1.6	11.1	9.3	0.09	6.6	8.6	8.9	4.6	90	192
C22	7.1	5.08	18	44.7	88	35.5	40.3	191	38.7	2.8	11.4	9.2	0.13	6.7	8.4	8.6	5	86	176

Table 6

Reference Patients	Sick	Healthy
Amlgm1	1	0
Amlgm2	1	0
Amlgm3	1	0
Amlgm4	1	0
Amlgm6	1	0
Amlgm7	1	0
Amlgm8	1	0
Amlgm9	1	0
Amlgm10	1	0
Amlgm11	1	0
Amlgm13	1	0
Amlgm14	1	0
Amlgm15	1	0
Amlgm16	1	0
Amlgm17	1	0
Amlgm20	1	0
Amlgm23	1	0
Amlgm24	1	0
Amigm25	1	0
Healthy Patients		-
aC1	0	. 1
aC2	0	1
aC3	0	1
aC4	0	1
aC5	0	1
aC7	0	1
aC8	0	1
C13	0	1
C14	0	1
C16	0	1
C18	0	1
C20	0	1
C15	0	1
C17	0	1
C21	0	1
C22	0	1

The second step involves the above described artificial neural network analytic model to select the matching disorder for the unknown and undiagnosed patient, and the associated PLS-Discriminant Analysis model, according to the above described process.

5

10

In the third and final step, predictions using the selected PLS-Discriminant Analysis model in step 2 are done in order to produce the diagnosis and the associated probability level. These are obtained by using one of the models (the selected model) developed in step 1. Again, the computation is carried out using Matlab® for Windows™ through use of the

program called "PLSPRED.M". 6 In this case, input is made to the equation -

ypred = plspred(x,b,P,Q,W,lv)

wherein

- ypred is the predicted value, i.e., the diagnosis; 5
  - x is the matrix of input vectors, comprised of data from undiagnosed patients;
  - b is the vector of regression coefficients between x and y, calculated in the first step;
- 10 • P is the matrix of x-loadings, calculated in the first step;
  - ullet Q is the matrix of y-loadings, calculated in the first step;
  - W is the matrix of weights, calculated in the first step;
- lv is the total number of latent variables, calculated in 15 the first step

function ypred = plspred(x,b,P,Q,W,lv)

- This program uses pls model parameters to predict the
- % responses to an new x block. Input/output sequence is: 20
  - ypred = plspred(x,b,P,Q,W,lv);
  - % Copyright
  - % Barry M. Wise
  - % 1991

25 [mx,nx] = size(x);

[mq, nq] = size(Q);

that = zeros(mx,lv);

ypred = zeros(mx,mq);

for i = 1:lv

30 that(:,i) = x\*W(:,i);

x = x - that(:,i)\*P(:,i)';

end

for i = 1:lv

ypred = ypred + b(1,i)\*that(:,i)\*Q(:,i)';

The computer screen indicates that the program is the copyright of Barry M. Wise, 1991.

33

end

After predictions have been obtained for the reference group in order to calculate **E**, the expected value in the equation below and the estimate of **s**, the standard deviation for the said reference group and for **X**, the predicted diagnosis, the t-value can be calculated as follows;

$$t = |X - E|/(s\sqrt{n})$$

wherein

- X is the value predicted by the PLS discriminant analysis
  - E is the expected value for the appropriate reference group;
  - s is the standard deviation associated with the reference group;
  - n is the number of observations.

The t-value obtained in this manner can be compared to tabulated values in a standard t-distribution and the associated probability interval determined.

20

25

30

35

15

5

As noted above, the invention is applicable to the analyses from essentially any method employed to determine These methods are well estabphysical medical disorders. lished and the literature abounds with descriptions of how they are conducted. It would be well within the skill of the art to employ the data from such analysis in the method of A commonly employed method of determining the invention. physical medical disorder is by immunoassay. In a typical immunoassay, the analyte is immunoreactive and its presence in a sample may be determined by virtue of its immunoreaction In a typical protein binding assay, with an assay reagent. the presence of analyte in a sample is determined by the specific binding reactivity of the analyte with an assay reagent where the reactivity is other than immunoreactivity. ples of this include enzyme-substrate recognition and the binding affinity of avidin for biotin. In the typical nucleic acid hybridization assay, the presence of analyte in a sample is determined by a hybridization reaction of the analyte with

10

15

20

25

30

35

WO 97/17891 PCT/US96/18254

an assay reagent. Analyte nucleic acid (usually present as double stranded DNA or RNA) is usually first converted to a single stranded form and immobilized onto a carrier (e.g., nitrocellulose paper). The analyte nucleic acid may alternatively be electrophoresed into a gel matrix. The immobilized analyte may then be hybridized (i.e., specifically bound) by a complementary sequence of nucleic acid.

The foregoing specific binding assays may be performed in a wide variety of assay formats. These assay formats fall within two broad categories. In the first category, the assay utilizes a chemiluminescent, bioluminescence, fluoroluminescence, and the like, or radio label conjugate which comprises the chemiluminescent, bioluminescent, fluoroluminescent, and the like, or the radio label moiety, attached to Specific binding material is a specific binding material. any material which will bind specifically by an immunoreaction, protein binding reaction, nucleic acid hybridization reaction, and any other reaction in which the material reacts specifically with a restricted class of biological, biochemi-In this category of assays, the cal or chemical species. conjugate participates in a specific binding reaction and the presence of analyte in the sample is proportional to the formation of one or more specific binding reaction products con-The assay is performed by allowing taining the conjugate. the requisite specific binding reactions to occur under suitable reaction conditions. The formation of specific binding reaction products containing the conjugate is determined by measuring the chemiluminescence, bioluminescence fluoroluminescence or energy emissions of such products containing the conjugate or by measuring the luminescent or radio emitting properties of unreacted or partially reacted conjugate not contained in such products.

This first category of assay formats is illustrated by sandwich assays, competitive assays, surface antigen assays, sequential saturation assays, competitive displacement assays and quenching assays.

In a sandwich format, the specific binding material to which the chemiluminescent, bioluminescent or fluoroluminescent moiety or radio label is attached, is capable of specifically binding with the analyte. The assay further utilizes a reactant which is capable of specifically binding with the analyte to form a reactant-analyte-labeled conjugate complex. The reactant may be attached to a solid phase, including without limitation, dip sticks, beads, tubes, paper or polymer sheets. In such cases, the presence of analyte in a sample will be proportional to the signal from the solid phase after the specific binding reactions are completed. Such assay formats are discussed further in U.S. Patent Nos. 4,652,533, 4,383,031, 4,380,580 and 4,226,993, which are incorporated herein by reference.

In a competitive format, the assay utilizes a reactant which is capable of specifically binding with the analyte to form an analyte-reactant complex and with the specific binding material, to which the labeled moiety is attached, to form a labeled conjugate-reactant complex. The reactant may be attached to a solid phase, or alternatively reaction products containing the reactant may be precipitated by use of a second antibody or by other known means. In this competitive format, the presence of analyte is "proportional," i.e., inversely proportional, to the signal from the solid phase or precipitate. A further discussion of this assay format may be found in the immediately above mentioned U.S. patents.

In another assay format, the analyte may occur on or be bound to a larger biological, biochemical or chemical species. This type of format is illustrated by a surface antigen assay. In this format, the specific binding material is capable of specifically binding with the analyte and the presence of analyte is proportional to the analyte-labeled conjugate complex formed as a reaction product. This is illustrated by attaching the label to an antibody which is specific to a surface antigen on a cell. The presence of the

cell surface antigen will be indicated by the signal from the cells after the completion of the reaction. The cells themselves may be used in conjunction with a filtration system to separate the analyte-label conjugate complex which is formed on the surface of the cell from unreacted labeled conjugate. This is discussed further in U.S. Patent No. 4,652,533.

36

PCT/US96/18254

WO 97/17891

5

10

15

20

25

30

35

Other assay formats include sequential saturation and competitive displacement, both of which utilize a labeled conjugate, where in both involve (1) the specific binding material, to which the moiety is attached, and (2) the analyte specifically bind with a reactant. In the case of sequential saturation, the analyte is reacted with the reactant first, followed by a reaction of the conjugate with remaining unreacted reactant. In the case of competitive displacement, the conjugate competitively displaces analyte which has already bound to the reactant.

In a quenching format, the assay utilizes a reactant which is capable of specifically binding with the analyte to form an analyte-reactant complex and with the specific binding material, to which the signal label is attached, to form a conjugate-reactant complex. A quenching moiety is attached to the reactant. When brought into close proximity to the label, the quenching moiety reduces or quenches the signal from the label. In this quenching format, the presence of analyte is proportional to the signal from the label. A further discussion of this format may be found in U.S. Patent Nos. 4,220,450 and 4,277,437, which are incorporated herein by reference.

In consideration of the above discussed assay formats, and in the formats to be discussed below, the order in which assay reagents are added and reacted may vary widely as is well known in the art. For example, in a sandwich assay, the reactant bound to a solid phase may be reacted with an analyte contained in a sample and after this reaction

the solid phase containing complexed analyte may be separated from the remaining sample. After this separation step, the label conjugate may be reacted with the complex on the solid phase. Alternatively, the solid phase, sample and label conjugate may be added together simultaneously and reacted prior to separation. As a still further but less preferred alternative, the analyte in the sample and the luminescent conjugate may be reacted prior to addition of the reactant on the solid phase. Similar variations in the mixing and reaction steps are possible for competitive assay formats as well as other formats known in the art.

In the second category of assay formats, the assay utilizes an unconjugated label compound. The presence of analyte in the sample is proportional to the formation of one or more specific binding reaction products which do not themselves contain the label. Instead, the label's signal is in proportion to the formation of such reaction products. example of this second category of assays, the assay utilizes a reactant capable of binding with the analyte to form an analyte-reactant complex which causes the signal from the la-This is illustrated by a simple enzyme-substrate assay in which the analyte is the substrate glucose and the reactant is the enzyme glucose oxidase. Formation of the enzymesubstrate complex triggers the luminescent compound. enzyme-substrate assay for glucose is disclosed in U.S. Patent Nos. 3,964,870 and 4,427,770, both of which are incorporated herein by reference. This enzyme-substrate assay is a specific binding assay in the sense that the substrate specifically binds to the active site of the enzyme in much the same way that an antigen binds to an antibody. In this assay, the enzyme specifically binds with the substrate which results in the production of peroxide which, in turn, causes the luminescent compound to luminesce.

35

10

15

20

25

30

Also included in the second category of assays are those assays in which the formation of the reaction products promotes or inhibits the signal by the label in a less direct

10

15

20

25

30

35

PCT/US96/18254 WO 97/17891

In this assay, a first reactant, which is cross reactive with the analyte, is attached to an enzyme such as glucose oxidase close to its active site. A second reactant which is specific for both the analyte and the immunoreactive material is added to the sample and the altered enzyme in the presence of the substrate (i.e., glucose). When the second reactant binds to the first reactant located near the active site on the enzyme, the second reactant blocks the active site in a way that the substrate cannot bind to the enzyme at the active site or the binding of the substrate at the active The second reactant blocksite is significantly decreased. ing the enzyme in this manner inhibits the enzyme from producing peroxide which, in turn, would have triggered the la-Analyte in the sample, however, will tie up the second reactant, thus preventing the second reactant from inhibiting the production of peroxide. The presence of analyte will be proportional to the signal from the compound.

The assays contained in the above two categories of assay formats may be heterogeneous or homogeneous. geneous assays, the reaction products, whose formation is proportional to the presence of analyte in the sample, are separated from other products of the reaction. Separation can be achieved by any means, including without limitation, separation of a liquid phase from a solid phase by filtration, microfiltration, double antibody precipitation, censize exclusion chromatography, removal of a trifugation, solid phase (e.g., a dip stick) from a sample solution or For example, in a sandwich assay the reacelectrophoresis. tant-analyte-label conjugate complex is separated from unreacted label conjugate. In a surface antigen assay, the analyte-label conjugate complex is separated form unreacted la-In a competitive assay, the reactant-label bel conjugate. conjugate complex is separated from unreacted label conjugate. In a sequential saturation assay and in a competitive displacement assay, the reactant-label conjugate complex is separated from unreacted label conjugate. Alternatively, in homogeneous assays the reaction products are not separated.

WO 97/17891

After the assay reagents have been allowed to react, the signal may be measured from the whole assay mixture whether such mixture is in solution, on a solid phase or distributed between various membrane layers of a dip stick or other solid support. The glucose assay using glucose oxidase and a label moiety illustrates a simple homogeneous assay in which separation is unnecessary. The quenching assay illustrates a more complex homogeneous assay in which separation is unnecessary. It is contemplated that either category of assay formats may give rise to either heterogeneous or homogeneous formats.

Other immunological antigen-antibody reactions methods suitable for generating the data for use in the method of the invention include, without limitation, turbidometry/nefelometry, electroimmunoassay (EIA), fluorescence polarized immunoassay (FPIA), electrophoresis, enzyme immunoassay (EIA), microparticle enzyme immunoassay (MEIA), and the like.

20

25

30

35

5

10

15

Wet and dry chemical analyses are an abundant source of data for the method of the invention. Analyses by large autoanalyzers can be used to generate such data. Chemical reactions that yield colored complexes, measured by spectrophotometry, for electrolytes, alkaline earth metal ions, and the like. Such can be effected also by the ion selective electrodes (ISE) method.

Measurements of enzymes such as alanine amino transferase, aspartate amino transferase, and the like; creatinin, urea, glucose, lactate, and drugs such as salicylate, paracetamol, and the like, are suitable data for the process of the invention.

Chromatographic analysis by HPLC, GC, GC/MS and liquid chromatography-mass spectrometry (LC/MS) can be used for data relative to serotonin metabolites, hemoglobin, catechol

amines, adrenaline, noradrenaline, dopamine, alcohols such as methanol, ethanol, n-propanol, and the like.

40

The data can be that generated by an amino acid sequenator relating to amino acid sequencing and DNA analysis, or by polymerase chain reaction (PCR) as used in forensic sciences, and coagulation techniques such as clotting analysis and chromogen substrate analysis.

5

10

#### CLAIMS:

5

10

15

20

25

30

1. A process for analyzing for a variety of physical medical disorders, which comprises:

- a) clinical testing by common analyses, patients free of certain physical medical disorders and patients that possess such one or more disorders, which analyses distinquish such disorders thereby obtaining the characterization of a collection of such disorders as numerical data;
- b) scaling the matrix of said numerical data in a computer;
  - c) configuring in a computer, an artificial neural network prescribed by the number of variables and the number of disorders, which artificial neural network is electronic data that possesses an input layer, one or more hidden layers and an output layer;
  - d) fitting the artificial neural network electronic data in a computer to the numerical data according to adjustable parameters of the artificial neural network;
  - e) whereby the artificial neural network is trained in the computer to automatically provide a analytic model;
  - f) withdrawing from a new patient that has not been diagnosed for such disorder, a sample or samples of a kind taken from said reference patients and subjecting said diagnostic sample to said clinical testing to obtain new numerical data;
  - g) automatically scaling in the computer the new data to the data derived from the reference patients;
  - feeding the scaled new data to the trained artificial neural network and the analytic model thereof, and
  - i) automatically obtaining a diagnosis of the new patient with respect to such disorders.
- 2. The process of claim 1 wherein the clinical testing involves diagnostic procedures that involves determining the nature and cause of a disease or injury through evaluation of patient history, examination, and review of laboratory data.

3. The process of claim 2 wherein the laboratory data includes the diagnostic analysis of bodily fluids, cells, soft and hard tissues, and organs.

- 4. The process of claim 3 wherein the data is from one or more of immunoassays, wet and dry chemical analysis, hematological analysis, cell size analysis, chromatographic analysis, spectrographic analysis, electrochemical analysis, pH analysis, ion selective potentiometry, NIR, ISS, SIMS, AES,
- ICP-MS, AAS, LC/MS, nuclear magnetic resonance spectrometry, nuclear magnetic imaging techniques, chemical elemental analysis, and gas analysis including gas chromatography, antibody analysis, blotting techniques, electrophoresis, amino acid sequencing, and DNA analysis.

5. The process of claim 1 wherein the obtained numerical data characterizes deviations from the data of healthy patients as well as data that correlates with the data of the reference group patients that possess the disorder.

20

35

6. The process of claim 1 wherein the obtained diagnosis provides the best matching reference group and thereby assigning the particular disorder for the physician.

7. The process of claim 1 wherein if no matching disorder is provided by the obtained diagnosis, any abnormality of the sample is flagged.

8. The process of claim 1 wherein autoscaling is done by analog or digital means and the translations are computed by computer.

9. The process of claim 8 wherein translating the data is effected on a computer through digital conversions.

10. The process of claim 1 wherein the artificial neural network is optimized by (i) determining the right number of hid-

97/17891 4.3

5

15

WO 97/17891 PCT/US96/18254

den layers, (ii) the right number of neurons in the hidden layers, and (iii) overall training of the network.

11. The process of claim 10 wherein (i) and (ii) are chosen by configuring the network by a pre-chosen number and training the network, calculating prediction error between real and model values, storing this data and choosing the next higher number for number (i) and (ii) above.

10 12. The process of claim 11 wherein the pre-chosen number is at least one hidden layer and at least two neurons.

13. The process of claim 12 wherein the new network is trained, prediction error is calculated and stored, and this process is repeated until a minimum prediction error is attained.

## INTERNATIONAL SEARCH REPORT

International application No. PCT/US96/18254

	SSIFICATION OF SUBJECT MATTER		
IPC(6) :	:A61B 5/00 -128/630		
	o International Patent Classification (IPC) or to both national classification and	IPC	
B. FIEL	DS SEARCHED		
Minimum de	ocumentation searched (classification system followed by classification symbols	)	
U.S. : 2	128/630, 653.1; 364/413.01, 413.02, 413.13; 395/21, 23, 24		
Documentat	ion searched other than minimum documentation to the extent that such document	s are included	in the fields searched
Electronic d	lata base consulted during the international search (name of data base and, whe	re practicable	, search terms used)
C. DOC	CUMENTS CONSIDERED TO BE RELEVANT		
Category*	Citation of document, with indication, where appropriate, of the relevant	passages	Relevant to claim No.
Υ	US 5,463,548 A (ASADA et al) 31 October 1995 line 44 to col. 5, line 17.	, col. 4,	1-13
Y	SHARPE et al., Artificial Neural Networks in Diag Thyroid Function from in Vitro Laboratory, Clinical Ch Vol. 39, No. 11, 1993, pages 2248-2253, especia 2249 col.1, lines 1-31.	emistry,	1-13
<b>Y</b>	WILDING et al., Application of backpropatgation networks to diagnosis of breast and ovarian cancer, Scientific Publishers Ireland Ltd., Cancer Letters 77, pages 145-153, especially page 146 col. 1, line 10 line 11.	Elsevier (1994),	1-13
Furth	ner documents are listed in the continuation of Box C. See patent far	nily annex.	
• Sp			ernational filing date or priority
	cument defining the general state of the art which is not considered principle or theory to be of particular relevance	inderlying the inv	ation but cited to understand the restion
	rlier document published on or after the international filing date "X" document of particles document published on or after the international filing date	cannot be conside	e claimed invention cannot be red to involve an inventive step
	cument which may throw doubts on priority claim(s) or which is when the document	is taken alone	
*O* do:	ccial reason (as specified)  "Y"  document of particle  considered to invo	ive an inventive ormore other suc	e claimed invention cannot be step when the document is h documents, such combination he art
*P* do	cument published prior to the international filing date but later than "&" document member of priority date claimed		
	actual completion of the international search  Date of mailing of the in	ternational ser	arch report
20 JANU	ARY 1997 17 MAR 19	97	
	mailing address of the ISA/US	این	C.
Box PCT Washington	n, D.C. 20231 FRANK JAWORSKI		
Facsimile N	lo. (703) 305-3590 / Telephone No. (703)	308-3961	<del></del>

Form PCT/ISA/210 (second sheet)(July 1992)\*